

### **Radiation Protection Aprons**

# Front Lead Apron with Buckle

Keeping you protected and comfortable during diagnostic imaging procedures with front exposure only





Closure A - Velcro flap between the shoulders



Closure B – Buckle with fully adjustable strap



Customised Just For YOU



Unique 3-Year Guarantee

#### **Product Features Overview**

- Radiation Protection Standards: IEC 61331.1-2014 & IEC 61331.3-2014, AU/NZS 4543.3-2000.
- With Standard Velcro Thyroid Collar and Pocket.
- Front Protection: one panel.
- · Buckle clip with fully adjustable strap for back support and even weight distribution.
- Upper Velcro strap to keep shoulder sleeves together.
- · Padded shoulders for more comfort.
- · Soft fabric material to allow for complete flexibility.

### **Custom Options**

- · Embroidery or Logo on pocket
- · Bib or Cosy Thyroid collar with Magnet or Velcro closure
- Extra pocket (e.g. dosimeter)
- Access to A-Track Inventory System
- · Available Fabrics:
  - Easy-Clean Flex-DC
  - Stain Resistant Spangly
  - Standard Nylon

<b>Protection Materials</b>	Lead Equivalence*	Protection	Size Options	
Lead	0.25		XS S M	
Light Lead Plus	0.35		Other sizes available on request.	
Supreme Lead Free	0.50	Front protection	Please refer to size chart on back page.	

<sup>\*</sup> Measurement in mm Pb.




## Sizing Chart for Front Lead Aprons

The following sizing chart caters for our most popular (standard) sizes. If you need help selecting the right size, please refer to the measuring guidelines below the sizing table.

Front Lead Aprons (in cm)						
Size	XS	S	М	L	XL	
Shoulder width	36	36	42	42	50	
Front width	50	55	60	75	75	
Standard length	100	100	100	100	100	

Shoulder width 36, 42 and 50 available on request. Other lengths available: 90, 110, 120. Other sizes available on request.



If your individual measurements do not fit the standard sizing charts below, please feel free to send them to info@deutschmedical.com.au and we will assist you in choosing the most suitable size for your radiation protection apron.

My Individual Measurements				
Distance shoulder to shoulder	cm			
Length from (highest point of) shoulders to below knee	cm			
Width at hip	cm			