

Radiation Protection Aprons

Front Lead Apron with Velcro

Keeping you protected and comfortable during diagnostic imaging procedures with front exposure only



Closure A – Velcro flap between the shoulders

Closure B – Elasticated velcro fastening provides individual fit

Velcro closure – side




Customised
Just For YOU



Unique 3-Year
Guarantee

Product Features Overview	Custom Options
<ul style="list-style-type: none"> • Radiation Protection Standards: IEC 61331.1-2014 & IEC 61331.3-2014, AU/NZS 4543.3-2000. • With Standard Velcro Thyroid Collar and Pocket. • Protection: one panel. • Adjustable Velcro flaps on hips for back support, even weight distribution and individual fit. • Upper Velcro strap to keep shoulder sleeves together. • Padded shoulders for more comfort. • Soft fabric material to allow for complete flexibility. 	<ul style="list-style-type: none"> • Embroidery or Logo on pocket • Bib or Cosy Thyroid collar with Magnet or Velcro closure • Extra pocket (e.g. dosimeter) • Access to A-Track Inventory System • Available Fabrics: <ul style="list-style-type: none"> - Easy-Clean Flex-DC - Stain Resistant Spangly - Standard Nylon

Protection Materials	Lead Equivalence*	Protection	Size Options
Lead	0.25	 <p>Front protection</p>	<div style="display: flex; justify-content: space-around;"> <div style="background-color: #0056b3; color: white; padding: 5px; border-radius: 10px;">XS</div> <div style="background-color: #0056b3; color: white; padding: 5px; border-radius: 10px;">S</div> <div style="background-color: #0056b3; color: white; padding: 5px; border-radius: 10px;">M</div> </div> <div style="display: flex; justify-content: space-around; margin-top: 5px;"> <div style="background-color: #0056b3; color: white; padding: 5px; border-radius: 10px;">L</div> <div style="background-color: #0056b3; color: white; padding: 5px; border-radius: 10px;">XL</div> <div style="background-color: #0056b3; color: white; padding: 5px; border-radius: 10px;">XXL</div> </div>
Light Lead Plus	0.35		Other sizes available on request.
Supreme Lead Free	0.50		Please refer to size chart on back page.


* Measurement in mm Pb.

Sizing Chart for Front Lead Aprons

The following sizing chart caters for our most popular (standard) sizes. If you need help selecting the right size, please refer to the measuring guidelines below the sizing table.

Front Lead Aprons (in cm)					
Size	XS	S	M	L	XL
Shoulder width	36	36	42	42	50
Front width	50	55	60	75	75
Standard length	100	100	100	100	100

Shoulder width 36, 42 and 50 available on request. Other lengths available: 90, 110, 120. Other sizes available on request.



If your individual measurements do not fit the standard sizing charts below, please feel free to send them to info@deutschmedical.com.au and we will assist you in choosing the most suitable size for your radiation protection apron.

My Individual Measurements	
Distance shoulder to shoulder	cm
Length from (highest point of) shoulders to below knee	cm
Width at hip	cm