

## Radiation Protection Aprons

# Full Wraparound Lead Apron

Keeping you protected and comfortable during diagnostic imaging procedures and surgeries



Product Features Overview	Custom Options
<ul style="list-style-type: none"> <li>• Radiation Protection Standards: IEC 61331.1-2014 &amp; IEC 61331.3-2014, AU/NZS 4543.3-2000.</li> <li>• Even transfer of weight between shoulders and hips to reduce fatigue.</li> <li>• With Standard Velcro Thyroid Collar and Pocket.</li> <li>• Front protection: two fully overlapping panels. Double-layer for shoulder area.</li> <li>• Back protection: one panel.</li> <li>• Adjustable Velcro for correct fit and padded shoulders for additional comfort.</li> <li>• Soft fabric material to allow for full flexibility.</li> </ul>	<ul style="list-style-type: none"> <li>• Embroidery or Logo on pocket</li> <li>• Bib or Cosy Thyroid collar with Magnet or Velcro closure</li> <li>• Extra pocket (e.g. dosimeter)</li> <li>• Access to A-Track Inventory System</li> <li>• Available Fabrics:               <ul style="list-style-type: none"> <li>- Easy-Clean Flex-DC</li> <li>- Stain Resistant Spangly</li> <li>- Standard Nylon</li> </ul> </li> </ul>

Protection Materials	Lead Equivalence*	Protection	Size Options
Lead	0.35/0.25    0.35/0.35	<p>Front, sides and back protection Full front overlapping</p>	<b>XS</b> <b>S</b> <b>M</b> <b>L</b> <b>XL</b> <b>XXL</b>
Light Lead Plus	0.50/0.25    0.50/0.35		Other sizes available on request.
Supreme Lead Free	0.70/0.25    0.70/0.35		Please refer to size chart on back page.


\* Measurement in mm Pb.

# Sizing Chart for Full Wraparound Lead Apron

The following sizing chart caters for our most popular (standard) sizes. If you need help selecting the right size, please refer to the measuring guidelines below the sizing table.

Full Wraparound Lead Apron (in cm)						
Size	XS	S	M	L	XL	XXL
	91	97	107	113	120	127
Chest circumference	80-85	86-91	92-97	98-107	108-113	114-120
Standard length	100	100	100	100	100	100

When measuring total length, allow for frontal overlap. For a better fit, measure at the widest point. Other lengths available: 90, 110, 120. Other sizes available on request.



If your individual measurements do not fit the standard sizing charts below, please feel free to send them to [info@deutschmedical.com.au](mailto:info@deutschmedical.com.au) and we will assist you in choosing the most suitable size for your radiation protection apron.

My Individual Measurements	
Circumference of chest at underarm	cm
Length from (highest point of) shoulders to below knee	cm