

Radiation Protection Aprons

Support Belt Lead Apron Skirt

Comfortable protection and lower back support during extended procedures



Customised
Just For YOU




Unique 3-Year
Guarantee

Product Features Overview

- Radiation Protection Standards: IEC 61331.1-2014 & IEC 61331.3-2014, AU/NZS 4543.3-2000.
- Supportive Wide Belt on internal back panel for lower back support.
- Front protection: two fully overlapping panels.
- Back protection: one panel.
- Adjustable Velcro closure at front for tight fit.
- Soft fabric material to allow for full flexibility.
- Inclusions: Wide Support Belt with Velcro.

Custom Options

- Embroidery of name/ department or logo on skirt
- Extra pocket (for dosimeter, or special)
- Access to A-Track Inventory Management Software
- Available Fabrics:
 - Easy-Clean Flex-DC
 - Stain Resistant Spangly
 - Standard Nylon

Protection Materials	Lead Equivalence*	Protection	Size Options
Lead	0.35/0.25 0.35/0.35	 <p>Front, sides and back protection Full front overlapping</p>	XS S M L LL XL XXL
Light Lead Plus	0.50/0.25 0.50/0.35		Other sizes available on request.
Supreme Lead Free	0.70/0.25 0.70/0.35		Please refer to size chart on back page.


* Measurement in mm Pb.

Sizing Chart for Vest & Skirt Lead Aprons

The following sizing chart caters for our most popular (standard) sizes. If you need help selecting the right size, please refer to the measuring guidelines below the sizing table.


Lead Apron Vest (in cm)						
Size	XS	S	M	L	XL	XXL
	91	97	107	113	120	127
Chest circumference	80-85	86-91	92-97	98-107	108-113	114-120
Length	50	55	60	65	65	65

When measuring total length allow for overlap. Other sizes available on request.



Lead Apron Skirt (in cm)							
Size	XS	S	M	L	LL	XL	XXL
Waist	65	70	75	81	87	93	103
Recommended waist measurements	63-67	68-72	73-78	79-84	85-90	91-97	98-108
Hips	90	95	100	106	112	118	128
Length	50	60	60	60	60	60	60

When measuring total length allow for overlap. Other lengths available: 70, 80. Any combination of size and length available. Other sizes available on request.



If your individual measurements do not fit the standard sizing charts below, please feel free to send them to info@deutschmedical.com.au and we will assist you in choosing the most suitable size for your radiation protection apron.

My Individual Vest Measurements	
Circumference of chest at underarm	cm
Length from (highest point of) shoulders to hip	cm

My Individual Skirt Measurements	
Waist circumference	cm
Hips circumference	cm
Length from waist to knee	cm