

## Radiation Protection Aprons

# Tailor-Fit Lead Apron Vest with Velcro

Keeping you fully protected and comfortable even during long procedures



Customised  
Just For YOU



Unique 3-Year  
Guarantee

Velcro closure

### Product Features Overview

- Radiation Protection Standards: IEC 61331.1-2014 & IEC 61331.3-2014, AU/NZS 4543.3-2000.
- Even transfer of weight between shoulders and hips to reduce fatigue.
- With Standard Velcro Thyroid Collar and Pocket.
- Front protection: two panels with 15cm overlap.
- Back protection: one panel.
- Adjustable Velcro at front and padded shoulders for comfort.
- Soft fabric material to allow for complete flexibility.

### Custom Options

- Embroidery or Logo on pocket
- Bib or Cosy Thyroid collar with Magnet or Velcro closure
- Extra pocket (e.g. dosimeter)
- Access to A-Track Inventory System
- Available Fabrics:
  - Easy-Clean Flex-DC
  - Stain Resistant Spangly
  - Standard Nylon

Protection Materials	Lead Equivalence*	Protection	Size Options
Lead	0.25/0.25	<p>Front, sides and back protection Partial front overlapping</p>	XS   S   M L   XL   XXL
Light Lead Plus	0.35/0.25   0.35/0.35		Other sizes available on request.
Supreme Lead Free	0.50/0.25   0.50/0.35		Please refer to size chart on back page.


\* Measurement in mm Pb.

# Sizing Chart for Vest & Skirt Lead Aprons

The following sizing chart caters for our most popular (standard) sizes. If you need help selecting the right size, please refer to the measuring guidelines below the sizing table.


Lead Apron Vest (in cm)						
Size	XS	S	M	L	XL	XXL
	91	97	107	113	120	127
Chest circumference	80-85	86-91	92-97	98-107	108-113	114-120
Length	50	55	60	65	65	65

When measuring total length allow for overlap. Other sizes available on request.



Lead Apron Skirt (in cm)							
Size	XS	S	M	L	LL	XL	XXL
Waist	65	70	75	81	87	93	103
Recommended waist measurements	63-67	68-72	73-78	79-84	85-90	91-97	98-108
Hips	90	95	100	106	112	118	128
Length	50	60	60	60	60	60	60

When measuring total length allow for overlap. Other lengths available: 70, 80. Any combination of size and length available. Other sizes available on request.



If your individual measurements do not fit the standard sizing charts below, please feel free to send them to [info@deutschmedical.com.au](mailto:info@deutschmedical.com.au) and we will assist you in choosing the most suitable size for your radiation protection apron.

My Individual Vest Measurements	
Circumference of chest at underarm	cm
Length from (highest point of) shoulders to hip	cm

My Individual Skirt Measurements	
Waist circumference	cm
Hips circumference	cm
Length from waist to knee	cm