

We offer a full range of preparation equipment, microscopes and accessories for material science and other applications

## Buehler Motopol 2000 Automatic Polishing machine



### Description

Buehler Motopol 2000 with platen for central pressure metallurgical sample preparation with automatic functions. The force for the machine is applied with compressed air and can be precisely controlled.

The base rotation speed is also set using the controller and offer contra and complimentary rotation. This machine is in great

### Features

- Power 220-240v 13amp
- Platen size 300mm
- Water inlet /outlet
- Compressed Air 4bar
- Weight 45KG
- Dimensions
- (L800mm x W500mm x H600mm)

**\*SAMPLE HOLDER IS AN OPTIONAL EXTRA AND THERE IS NO SPLASH RING**



[www.spectrographic.co.uk](http://www.spectrographic.co.uk)

To speak to a specialist for expert advice, call 01943 879001