

We offer a full range of preparation equipment, microscopes and accessories for material science and other applications

PetroTrim



Description

- The PetroTrim is an easy to use precision saw machine, it is used for the sectioning of samples.
- This precision saw offers a number of advantages over other saws. Particularly that the linear speed is of a variable system and the feed rate is monitored which allows the operator to control the rate at which the samples move through the blade. This evades many of the problems that occur when using gravity-fed systems. E.g. slides shattering if the entire cutting force is exerted on a smaller area of the specimen at the end of the cut.
- Stock Code: 910-030

Features

- Easy to use
- Simple interface
- Accurate cutting
- Disc Size- 150mm
- Arbor Size- 16mm
- Vice size- N/A
- Coolant connection size- 6mm pushfit
- Drain size- 20mm ID 24mm OD
- Coolant tank size- N/A
- Power- 240V
- Dimensions- 670x440x370 mm
- Weight- 20kg



www.spectrographic.co.uk

To speak to a specialist for expert advice, call 01943 879001