

We offer a full range of preparation equipment, microscopes and accessories for material science and other applications

EVO 400



Description

- The EVO 400 is a cutting machine that brings a cutting solution for large samples (up to 200mm high and over 200kg).
- This new machine is great for workshops and combines sturdiness with ease of use.
- Its programmable interface allows for repetitive cutting and has been designed with ergonomics in mind with a motorized hood and touchscreen control panel for easy and efficient operation.

Features

- Colour touch screen interface
- Password protection can be implemented on the machine
- Up to 100 password-protected programs
- LED internal lighting
- Integrated spray nozzle for cleaning
- Optional external 200L re-circulation tank
- Interlock switch keeps hood locked during cutting
- Rotation speed- 1000-3000RPM
- Cut off wheel size- Up to 400mm
- Power Supply- 380V, 3 Phases, 50Hz
- Motor Power- 7.5Kw
- Dimensions- 1140x1600x1910mm
- Weight 750Kg



www.spectrographic.co.uk

To speak to a specialist for expert advice, call 01943 879001