

Compatibility Chart

Description	Cat. No.	AHN Pipet4u Pro / Pipet4u Revolution single and multi channel pipettes				Gilson Pipetman / Pipetman G single and multi channel pipettes				Eppendorf Research Plus single and multi channel pipettes				Brand Transferpette, Transferpette S single and multi channel pipettes				Thermo Finnpipette F1 / Novus manual single and multi channel and electronic single and multi channel pipettes				Sartorius Praline Plus / Picuss manual single channel and electronic single and multi channel pipettes				HTL Labmate Pro (same as Labnet BioPette) manual single and multi channel pipettes											
		0.2-2 µL	0.5-10 µL	2-20 µL	5-50 µL	10-100 µL	20-200 µL	100-1000 µL	1000-5000 µL	0.5-10 µL	10-100 µL	20-200 µL	100-1000 µL	1000-10000 µL	0.2-2 µL	0.5-10 µL	2-20 µL	10-100 µL	20-200 µL	100-1000 µL	1000-10000 µL	0.5-10 µL	10-100 µL	20-200 µL	100-1000 µL	1000-10000 µL	0.5-10 µL	10-100 µL	20-200 µL	100-1000 µL	1000-10000 µL	0.5-10 µL	10-100 µL	20-200 µL	100-1000 µL	1000-10000 µL	
AHN myTip® Pipette Tips																																					
AHN myTip® 10 µL, clear	1-002-81-0	■	■																																		
AHN myTip® 10 µL, clear, sterile	1-003-C5-0	■	■																																		
AHN myTip® 10/20 µL XL, clear	1-062-81-0																																				
AHN myTip® 10/20 µL XL, clear, sterile	1-062-C5-0	■	■																																		
AHN myTip® 10/20 µL XL, clear, sterile	1-063-C5-0	■	■																																		
AHN myTip® 200 µL, clear	1-121-86-0																																				
AHN myTip® 200 µL, clear, sterile	1-121-C5-0																																				
AHN myTip® 200 µL, clear, sterile	1-122-C5-0																																				
AHN myTip® 200 µL, yellow	1-121-86-1																																				
AHN myTip® 200 µL, yellow, sterile	1-121-C5-1																																				
AHN myTip® 200 µL, yellow, sterile	1-122-C5-1																																				
AHN myTip® 200 µL, Gelloader, clear	1-130-70-0																																				
AHN myTip® 200 µL, Gelloader, clear	1-130-C4-0																																				
AHN myTip® 200 µL, Gelloader, clear, sterile	1-131-C4-0																																				
AHN myTip® 300 µL, clear	1-127-80-0																																				
AHN myTip® 300 µL, clear, sterile	1-127-C5-0																																				
AHN myTip® 300 µL, clear, sterile	1-128-C5-0																																				
AHN myTip® 1000 µL, short, clear blue	1-209-80-2																																				
AHN myTip® 1000 µL, clear	1-202-80-0																																				
AHN myTip® 1000 µL, clear	1-202-C4-0																																				
AHN myTip® 1000 µL, clear, sterile	1-203-C4-0																																				
AHN myTip® 1000 µL, blue	1-202-80-2																																				
AHN myTip® 1000 µL, blue, sterile	1-202-C4-2																																				
AHN myTip® 1000 µL, blue, sterile	1-203-C4-2																																				
AHN myTip® 1250 µL, clear	1-204-80-0																																				
AHN myTip® 1250 µL, clear, sterile	1-204-C4-0																																				
AHN myTip® 1250 µL, clear, sterile	1-205-C4-0																																				
AHN myTip® Mat Macro Pipette Tips																																					
AHN myTip® Mat 5 mL, clear	1-218-61-0																																				
AHN myTip® Mat 5 mL, clear, sterile	1-222-61-0																																				
AHN myTip® Mat 5 mL, clear	1-210-84-0																																				
AHN myTip® Mat 5.5 mL, clear	1-212-84-0																																				
AHN myTip® Mat 5.5 mL, clear, sterile	1-216-84-0																																				
AHN myTip® Mat 10 mL, clear	1-213-50-0																																				
AHN myTip® Mat 10 mL, clear	1-220-50-0																																				
AHN myTip® Mat 10 mL, clear, sterile	1-224-50-0																																				

■ compatible x not compatible

The compatibility chart above is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes. Specifications like inaccuracy and imprecision values have not been tested. For additional compatibility charts, please contact info@cappahn.com.



+49 (0) 3631 / 65242-0
 +49 (0) 3631 / 65242-90
 Uhlheber Weg 14 | 99734 Nordhausen | Germany
info@cappahn.com
www.dhn-bio.de



