

## Sling Fabric Selection Chart

### Disposable



- Single client use
- Non-woven
- Stitch-bond fabric
- Liquid repellent
- Light soiling can be wiped
- Semi-air permeable
- Should not be laundered

### PolySlip



- Closed weave polyester knit
- Comfortable
- Durable
- High strength
- Breathable
- Washable

### Quilted Disposable



- Single client use
- Padded
- Non-woven
- Increased comfort
- Reduces bunching
- Liquid repellent
- Light soiling can be wiped
- Should not be laundered

### Slip Fit



- Nylon ripstop dual layer material
- Easily fitted and removed
- High-strength
- Durable
- Non-breathable
- Washable

### Mesh



- Open-weave polyester
- High strength
- Air permeable
- Can be used for bathing
- Durable
- Washable

### Spacer



- Knitted polyester with an internal air gap
- Multi-directional stretching
- Used for bathing
- Air permeable
- Pressure redistributing
- Can be left behind patient
- Washable

### Quilted Mesh



- Mesh fabric with quilted foam and polyester
- Air permeable
- Increased comfort
- Reduces bunching
- Washable

**Please note:** not all slings are available in all fabric options. Please refer to the website for available options.