

TiLobe® Multi-Unit Abutment System

The TiLobe® Multi-Unit Abutment is designed for screw-retained restorations of partially edentulous and fully edentulous arches, such as the Keystone Arch treatment concept. These abutments are available in both straight (non-engaging, one-piece) and angled (17° and 30°) alternatives with a selection of cuff heights to match soft tissue thickness.

Material: **Titanium Alloy**


Recommended Abutment/
Abutment Screw Torque Value: **30 Ncm**

Recommended Prosthetic Screw
Torque Value: **20 Ncm**

Driver Required: **Quad Driver**


Indications For Use:

- Multiple-unit screw retained restorations
- Edentulous or partially edentulous arches
- Simplified prosthetic access for diverging implant angles
- Transition from table of implant through the soft tissue



Prosthetic Components

Description	Catalog #
Multi Unit Abutment Healing Cap	45330K
Multi Unit Abutment Impression Post, Open Tray (with screw)	45327K
Multi Unit Abutment Impression Post, Closed Tray	45328K
Multi Unit Abutment Analog	45329K
Multi Unit Abutment Ti Temporary Sleeve (with prosthetic screw)	45331K
Multi Unit Abutment Plastic Sleeve (with prosthetic screw)	45332K
Replacement Multi Unit Abutment Prosthetic Screw	45333K
Replacement Abutment Screw for Angled Multi Unit Abutment	45334K



Ø3.5/3.8mm Internal TiLobe® Implants

Description	Cuff	Catalog #
Straight	1.0 mm	45312K
Straight	2.0 mm	45313K
Straight	3.0 mm	45314K
Angled 17°	3.0 mm	45321K
Angled 17°	4.0 mm	45322K
Angled 30°	3.0 mm	45936K
Angled 30°	4.0 mm	45937K
Angled 30°	5.0 mm	45938K

Ø4.1/4.5mm Internal TiLobe® Implants

Description	Cuff	Catalog #
Straight	1.0 mm	45315K
Straight	2.0 mm	45316K
Straight	3.0 mm	45317K
Straight	4.0 mm	45318K
Angled 17°	3.0 mm	45323K
Angled 17°	4.0 mm	45324K
Angled 30°	3.0 mm	45965K
Angled 30°	4.0 mm	45325K
Angled 30°	5.0 mm	45326K

Ø5.0/5.5/6.0/6.5mm Internal TiLobe® Implants

Description	Cuff	Catalog #
Straight	1.0 mm	45319K
Straight	2.0 mm	45320K
Straight	3.0 mm	45939K
Straight	4.0 mm	45940K
Angled 17°	3.0 mm	45941K
Angled 17°	4.0 mm	45942K
Angled 30°	3.0 mm	45966K
Angled 30°	4.0 mm	45967K

Quad Drivers

Description	Catalog #
Quad Driver, Latch, X-Long, 35 mm	45294K
Quad Driver, Assy, Swivel, X-Long, 35 mm	45293K

