

TITANIUM BASES ROUND

The Ti-Base for individual abutments or screw retained crown & bridge serves as a connector between the implant and its final restoration, enabling CAD/CAM customized solutions, with high precision and excellent esthetic results.

Paltop Ti-Bases of variable gingival heights provide solutions for cases of different gingival thickness and clinical procedures. The titanium base has Paltop's unique concave design enabling the soft tissue to form a dense layer surrounding the transgingival interconnection.

Short gingival height will be advantageous when the final restoration is close to the implant level, while long gingival height will be useful in cases when a high emergence profile is required (e.g. Sub-crestal placement).

Our new Round Ti-bases provide an excellent prosthetic CAD/CAM solution.

Ti-bases Round are available in a number of design options:



1. Titanium Bases Round

(recommended for posterior teeth).
4mm length are compatible with Sirona
and can be used with Sirona Scanbody.



2. Titanium Bases Round Anterior

(recommended for anterior teeth)

TITANIUM BASES ROUND

NARROW PLATFORM

Engage	LENGTH = 4 MM			LENGTH = 6 MM		
	0.5mm	2mm	3mm	0.5mm	2mm	3mm
	40-70215	40-70295	40-70217	40-70219	40-70297	40-70221
Non Engage	LENGTH = 4 MM			LENGTH = 6 MM		
	0.5mm	2mm	3mm	0.5mm	2mm	3mm
	40-70216	40-70296	40-70218	40-70220	40-70298	40-70222

STANDARD PLATFORM

Engage	LENGTH = 4 MM			LENGTH = 6 MM		
	0.5mm	2mm	3mm	0.5mm	2mm	3mm
	40-70207	40-70291	40-70209	40-70211	40-70293	40-70213
Non Engage	LENGTH = 4 MM			LENGTH = 6 MM		
	0.5mm	2mm	3mm	0.5mm	2mm	3mm
	40-70208	40-70292	40-70210	40-70212	40-70294	40-70214

WIDE PLATFORM

Engage	LENGTH = 4 MM	LENGTH = 6 MM	Non-Engage	LENGTH = 4 MM	LENGTH = 6 MM
	0.5mm	0.5mm		0.5mm	0.5mm
	40-70223	40-70225		40-70224	40-70226

CONICAL PLATFORM

Engage	LENGTH = 4 MM			LENGTH = 6 MM			Non Engage	LENGTH = 4 MM			LENGTH = 6 MM		
	0.5mm	2mm	3mm	0.5mm	2mm	3mm		0.5mm	2mm	3mm	0.5mm	2mm	3mm
	40-72046	40-72084	40-72048	40-72050	40-72086	40-72052		40-72047	40-72085	40-72049	40-72051	40-72087	40-72053

TITANIUM BASES ROUND ANTERIOR

NARROW PLATFORM STANDARD PLATFORM CONICAL PLATFORM

Engage	LENGTH = 4 MM			Engage	LENGTH = 4 MM			Engage	LENGTH = 4 MM		
	0.5mm	2mm	3mm		0.5mm	2mm	3mm		0.5mm	2mm	3mm
	40-70302	40-70303	40-70304		40-70299	40-70300	40-70301		40-72081	40-72082	40-72083