

## ALIGNER & RETAINER CUTTING, TRIMMING & FINISHING VIDEO

### **A Step-by-Step Video Demonstration**

We regularly interview our lab tech customers about the best practices and techniques for fabricating top-tier Zendura appliances and are pleased to share the video below with you. With Dr. Alejandro Ramírez conducting the demonstration, he uses two rotary cutting bits and one polishing head on an NSK rotary tool to cut the retainer from the Zendura sheet, trim away excess material, and then polish the edges of the final cutline.



We want to emphasize that this is one of several techniques for cutting, trimming and finishing Zendura appliances.

If you have special technique or procedure you think is worth sharing, or have questions, don't hesitate to contact our Customer Service department at 1-650-566-0800. Or email us at [info@ZenduraDental.com](mailto:info@ZenduraDental.com).

For an illustrated step-by-stem guide on how to cut, trim and finish Zendura appliances, [click here](#).