


# Zendura<sup>®</sup> AT

Attachment Template Material

Zendura AT™ is the first and only thermoformable plastic sheet material engineered from the start specifically for fabricating clear aligner therapy (CAT) attachment templates.

Now clinicians can form attachment templates that capture the attachment details designed in software more exactly than possible using other materials. Zendura AT also allows the clinician to more precisely place attachments, as well more easily remove the template once the attachments are placed.

Capture attachment detail more exactly.



Significantly more precise attachment placement.

A premium polyurethane material offering the perfect blend of rigidity, elasticity and clarity, Zendura AT is the latest material from Bay Materials, the Silicon Valley polymer R&D company recognized as the leading developer and supplier of innovative CAT materials. Zendura AT supports many of the most popular dental composite materials.

**BETTER TEETH  
ADAPTATION**

**BETTER ATTACHMENT  
DETAIL**

**MORE PRECISE  
PLACEMENT**

**NOTEABLY  
THINNER**

**BETTER  
CLARITY**

**BayMaterials<sup>®</sup>**

[www.zenduradental.com](http://www.zenduradental.com)

**+1.650.566.0800**



# Zendura<sup>®</sup> AT


Attachment Template Material

## The New Benchmark for Accuracy and Precision

- Better teeth adaptation
- Better attachment detail
- More precise attachment placement
- Flexible enough to allow easy removal of template once attachments are placed
- Noticeably thinner
- Better clarity
- Easy fabrication
- Supports many of the most popular dental composite materials

## The Difference Is Clear to See


Aluminum test model




**Zendura AT compared with another commercially available plastic sheet material suitable for making attachment templates.** Zendura AT conforms significantly more intimately to the 0.022" cut grooves on the top of the aluminum test model, following every detail of this challenging thermoforming adaptability test!




0.022" grooves cut into aluminum model




Competitive thermoformed sheet



Zendura AT



Zendura AT forms sharp details – modeling attachment and placement precisely



Other materials can form imprecisely ... with less precise location

### Order information

SKU# 9206 - Zendura AT 0.38mm (0.015") x 125mm

### For further information contact:

Bay Materials, Inc.  
 48450 Lakeview Blvd., Fremont, CA 94538  
 +1.650.566.0800 | Email: [info@zenduradental.com](mailto:info@zenduradental.com)  
[www.zenduradental.com](http://www.zenduradental.com)