

# Extensions - Partial Foot Toe Filler

Information Provided By Kevin Orthopedic Institute

## Toe Filler

Partial foot toe filler of Myolite, EVA and Plastazote material

### FUNCTION:

- Stabilizes anatomical segments adjacent to amputations

### CLINICAL INDICATION:

- Amputation

Toe filler is a wedge buttress of EVA, Myolite and Plastazote that sits on the superior forefoot area of a device's top cover. It conforms to and occupies the negative space from a patient's amputation(s).

**Note:** Submission of patient's shoe with impression is required for optimal fitting.



Toe Filler




5<sup>th</sup>

4<sup>th</sup>

3<sup>rd</sup>

2<sup>nd</sup>

1<sup>st</sup>



Superior view

**Notes:** All illustrations and diagrams are of right foot  
Colors on illustrations are for visual purposes and will vary on final product