## Offloading Pads & Cushions - Rearfoot

Information Provided By KevinRoot Medical



## Extrinsic rounded

Extrinsic rounded of frame rearfoot post

## **FUNCTIONS:**

- · Increased stability and strobel lasted shoes
- · Stabilizes rearfoot

## **CLINICAL INDICATIONS:**

- Overpronation
- Subtalar joint eversion/inversion
- · Shoe gear with strobel last

An extrinsic post that has been beveled along its plantar edges medially, proximally, and laterally. This functions to help fit into a variety of shoes including shoes with both slip-lasted and board-lasted construction (a rigid shoe platform made of plastic or cardboard).





Medial view







Varus Extrinsic shape for comparison

**Notes:** All illustrations and diagrams are of right foot Colors on illustrations are for visual purposes and will vary on final product ORF = Orthotic Reaction Force

call us: 1-800-496-0987

fax: 1-866-919-9268

