



Detached Carbon Foot Plate to Toes

Detached Carbon Foot Plate to Toes

Rigid Semi-rigid



Independent of orthotic device, carbon material extending the full length of foot model: heel to toes

FUNCTIONS:

- Immobilizes forefoot
- Reduces forefoot stress

CLINICAL INDICATIONS:

- Turf toe
- Hallux limitus
- Hallux rigidus
- Metatarsalgia
- Plantar plate injuries
- Arthritis

The detached carbon foot plate to toes is a 1-2mm thick solid piece of carbon material that extends from the heel and terminates just proximal to the end of the toes. This modification sits under the orthotic device between the shoe and orthotic.



Superior view
(transparent)



Plate fabricated to fit
patient's foot



Oblique lateral view



Lateral view

Notes: All illustrations and diagrams are of right foot
Colors on illustrations are for visual purposes and will vary on final product