

www.kevinrootmedical.com  
hello@kevinrootmedical.com  
Tel: 1 800 496 0987  
Fax: 1 866 919 9268

**Order Form**

Date MM...../DD...../YYYY.....

**Clinician**

Account Name/Number .....

Location .....

PO Number .....

Clinician .....

Contact me before processing

**Patient**

First Name .....

Last Name .....

Gender: M F DOB MM...../DD...../YYYY.....

Height ..... Weight ..... Shoe size .....

Dx .....

Ship to Patient

Street Address .....

City .....

State ..... Zip .....

**G140 Tall Gauntlet**



- 12" height
- 4mm polypropylene frame
- Rearfoot intrinsic post
- Leather top cover
- Leather liner
- Velcro and lace closures
- Device to mets

Suggested L Codes:  
L1940, L2820, L2330

**Impression**

- Plaster
- STS Casting Socks
- 3D Foot Scanner
- Fiberglass Casting Tape

**Fitting**

Fit AFO to submitted:

- Tracing
- Shoes

**Side**

Use separate Rx for each side

\*G800 exempt BIL standard

- Right
- Left

**Impression preparation**


**Ankle** As is Correct to 90°

**Rearfoot** As is Correct to 90°

**Forefoot** As is Balance FF to RF

**Measurements (optional)**

- Height:**
- Fibula Head
  - Proximal Trim
- Circumference:**
- Forefoot



- Circumference:**
- Proximal Trim
  - Above Ankle
  - Ankle
  - Mid-Foot

**Special instructions**

**Frame trim lines and options**

**Medial flange height**

**Lateral flange height**

**Closure options**

**Color options**

**Medial tibial trimlines**

**Lateral tibial trimlines**

**Frame foot plate**

**Heel cut out** Most stability Fits in shoe easier

**Pads**

- Met balance pad
- Met punch
- Metatarsal pad 2-4
- Met Bar 1-5
- Dancer's Pad
- Scaphoid Pad
- Cuboid Offload
- Heel pad
- Heel Spur Pad
- Plastazote liner (Lining under leather)

**Forefoot extensions**

**Sulcus extension**

- Sulcus Extension **Myolite (mm)**
- Morton's Ext.
- Rev. Morton's Ext.

**Toe extension**

- Morton's Ext. **Myolite (mm)**
- Rev. Morton's Ext.
- Toe Extension **Myolite (mm)**
- Toe Filler **Myolite (mm)**

**Frame modifications**

- 1<sup>st</sup> Ray Cut Out
- 5<sup>th</sup> Ray Cut Out
- Navicular B. Out (mm) (mm)
- 5<sup>th</sup> Button Out (mm) (mm)
- Heel Lift (mm) (mm)
- Heel Lift Tapered (mm) (mm)