

www.kevinrootmedical.com
hello@kevinrootmedical.com
Tel: 1 800 496 0987
Fax: 1 866 919 9268

Clinician

Account Name/Number
Location
PO Number
Clinician
Contact me before processing

G250 Overlapping Articulating Gauntlet



- 9" height, (5" and 12" optional)
- 4mm polypropylene frame
- Rearfoot intrinsic post
- Leather top cover, leather liner
- 2 anterior/1 posterior straps
- Exterior struts
- Device to mets

Suggested L Codes:
L1970, L2275, L2820, L2330

Impression

- Plaster
- STS Casting Socks
- 3D Foot Scanner
- Fiberglass Casting Tape

Fitting

- Fit AFO to submitted:
- Tracing
 - Shoes

Side

Use separate Rx for each side
*G800 exempt BIL standard



Impression preparation


Ankle As is Correct to 90°

Rearfoot As is Correct to 90°

Forefoot As is Balance FF to RF

Measurements (optional)

Height:
Fibula Head
Proximal Trim
Circumference:
Forefoot



Circumference:
Proximal Trim
Above Ankle
Ankle
Mid-Foot

Special instructions

Patient

First Name
Last Name
Gender: M F DOB MM...../DD...../YYYY.....
Height Weight Shoe size
Dx
Ship to Patient
Street Address
City
State Zip

Frame trim lines and options

Height	G109 - Standard (mid fibula)	G105 - Short (above malleoli)	G112 - Tall (to fibula head)
Medial flange height			
Lateral flange height			
Color options			

Pads

- Met balance pad
- Met punch
- Metatarsal pad 2-4
- Met Bar 1-5
- Dancer's Pad
- Scaphoid Pad
- Cuboid Offload
- Heel pad
- Heel Spur Pad
- Plastazote liner
(Lining under leather)

METATARSAL HEAD
(select all that apply)

Forefoot extensions

Sulcus extension		Sulcus Extension	Myolite (mm)	Myolite (mm)
		Morton's Ext.		
		Rev. Morton's Ext.		
Toe extension		Morton's Ext.		
		Rev. Morton's Ext.		
		Toe Extension	Myolite (mm)	Myolite (mm)
		Toe Filler		

Frame modifications

	1st Ray Cut Out		
	5th Ray Cut Out		
	Navicular B. Out	(mm)	(mm)
	5th Button Out	(mm)	(mm)
	Heel Lift	(mm)	(mm)
	Heel Lift Tapered	(mm)	(mm)