

www.kevinrootmedical.com  
hello@kevinrootmedical.com  
Tel: 1 800 496 0987  
Fax: 1 866 919 9268

**Clinician**

Account Name/Number .....  
Location .....  
PO Number .....  
Clinician .....  
Contact me before processing

**Patient**

First Name .....  
Last Name .....  
Gender: M F DOB MM...../DD...../YYYY.....  
Height ..... Weight ..... Shoe size .....  
Dx .....  
Ship to Patient .....  
Street Address .....  
City .....  
State ..... Zip .....

**G300 Tamarack Dorsi Assist Gauntlet**



- 9" height
- 4mm polypropylene frame
- Rearfoot intrinsic post
- Tamarack dorsi-assist joints
- Leather top cover, leather liner
- Optional closures
- Device to mets

Suggested L Codes:  
L1970, L2820, L2330, L2210

**Frame trim lines and options**

**Medial flange height** **Lateral flange height**

**Color options**

**Medial tibial trimlines** **Lateral tibial trimlines**

**Heel cut out** **Tamarack Dorsi-Assist Durometer options:** **Closure options**

**Impression**

- Plaster
- STS Casting Socks
- 3D Foot Scanner
- Fiberglass Casting Tape

**Fitting**

- Fit AFO to submitted:
- Tracing
  - Shoes

**Side**

Use separate Rx for each side  
\*G800 exempt BIL standard



**Impression preparation**

**Ankle** **Rearfoot** **Forefoot**

**Pads**

- Met balance pad
- Met punch
- Metatarsal pad 2-4
- Met Bar 1-5
- Dancer's Pad
- Scaphoid Pad
- Cuboid Offload
- Heel pad
- Heel Spur Pad
- Plastazote liner (Lining under leather)

**METATARSAL HEAD**  
(select all that apply)

**Forefoot extensions**

**Sulcus extension**

- Sulcus Extension **Myolite (mm)**
- Morton's Ext. **Myolite (mm)**
- Rev. Morton's Ext. **Myolite (mm)**

**Toe extension**

- Morton's Ext. **Myolite (mm)**
- Rev. Morton's Ext. **Myolite (mm)**
- Toe Extension **Myolite (mm)**
- Toe Filler **Myolite (mm)**

**Measurements (optional)**

**Height:** Fibula Head, Proximal Trim  
**Circumference:** Forefoot

**Circumference:** Proximal Trim, Above Ankle, Ankle, Mid-Foot

**Special instructions**

**Frame modifications**

- 1<sup>st</sup> Ray Cut Out (mm)
- 5<sup>th</sup> Ray Cut Out (mm)
- Navicular B. Out (mm)
- 5<sup>th</sup> Button Out (mm)
- Heel Lift (mm)
- Heel Lift Tapered (mm)