

1 CLINICIAN REQUIRED

Account Location

PO Number

Clinician

Clinician Email

A300 Tamarack Dorsi-Assist AFO Order Form

Date: / /
MM DD YYYY

Rush order due date:

..... / /
MM DD YYYY

Contact me to review Order Form



Call us: 1-877-767-3338

Fax: 1-866-919-9268

www.kevinorthopedic.com

hello@kevinorthopedic.com

2 PATIENT REQUIRED

Patient's Email

First Name

Last Name

DOB / / Sex: M F Weight Height



A300 Tamarack Dorsi-Assist AFO

- 4mm Polypropylene frame
- Rearfoot intrinsic post
- Tamarack dorsi joints
- 3mm 30 shore A EVA top cover
- Posterior bar
- Upper lining
- Mid fibula height
- Device to mets

3 SIDE

Use separate Rx for each side

Right

Left

IMPRESSION

Plaster

STS Casting Socks

3D Foot Scanner

Fiberglass Casting Tape

FITTING

Shoe size:

Fit AFO to submitted:

Tracing

Shoes


4 MODEL HEIGHT & COVER LENGTH REQUIRED

Height:

A109 - Standard (mid fibula)

A105 - Short (above malleoli)

A112 - Tall (to fibula head)



Length of 3mm EVA Cover:

To Mets

To Sulcus

To Toes

A MEASUREMENTS OPTIONAL

Height:

Fibula Head

Proximal Trim

Circumference:

Forefoot

Circumference:

Proximal Trim

Above Ankle

Ankle

Mid-Foot

B PRONATION CORRECTIONS OPTIONAL

	Varus Forefoot Posts	Balance FF to RF	8° <input type="checkbox"/> 6° <input type="checkbox"/> 4° <input type="checkbox"/> 2° <input type="checkbox"/> 0° <input type="checkbox"/> 90° <input type="checkbox"/>	
		VR Extrinsic (Bar)	6° <input type="checkbox"/> 4° <input type="checkbox"/> 2° <input type="checkbox"/> 0° <input type="checkbox"/>	
		VR Met Wedge	6° <input type="checkbox"/> 4° <input type="checkbox"/> 2° <input type="checkbox"/>	
	Midfoot	Raise Med Arch	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Scaphoid Pad	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		VR Cuboid Pad	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Medial Flange	High <input type="checkbox"/> Medium <input checked="" type="checkbox"/> Low <input type="checkbox"/>	
	Varus Rearfoot Posts	VR Intrinsic	6° <input type="checkbox"/> 4° <input type="checkbox"/> 2° <input type="checkbox"/> 0° <input checked="" type="checkbox"/>	
		VR Extrinsic	6° <input type="checkbox"/> 4° <input type="checkbox"/> 2° <input type="checkbox"/> 0° <input type="checkbox"/>	
		Heel Skive	6 <input type="checkbox"/> 4 <input type="checkbox"/> 2 <input type="checkbox"/> (mm)	

C SUPINATION CORRECTIONS OPTIONAL

	Valgus Forefoot Posts	Balance FF to RF	8° <input type="checkbox"/> 6° <input type="checkbox"/> 4° <input type="checkbox"/> 2° <input type="checkbox"/> 0° <input type="checkbox"/> 90° <input type="checkbox"/>	
		VG Extrinsic (Bar)	6° <input type="checkbox"/> 4° <input type="checkbox"/> 2° <input type="checkbox"/> 0° <input type="checkbox"/>	
		VG Met Wedge	6° <input type="checkbox"/> 4° <input type="checkbox"/> 2° <input type="checkbox"/>	
	Midfoot	Lower Med Arch	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Lateral Flange	High <input type="checkbox"/> Medium <input type="checkbox"/> Low <input type="checkbox"/>	
		VG Frame Filler	6° <input type="checkbox"/> 4° <input type="checkbox"/> 2° <input type="checkbox"/> 0° <input type="checkbox"/>	
		VG Cuboid Pad	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
	Valgus Rearfoot Posts	VG Intrinsic	6° <input type="checkbox"/> 4° <input type="checkbox"/> 2° <input type="checkbox"/> 0° <input type="checkbox"/>	
		VG Extrinsic	6° <input type="checkbox"/> 4° <input type="checkbox"/> 2° <input type="checkbox"/> 0° <input type="checkbox"/>	

D EXTENSIONS OPTIONAL

	Toe Length	Toe Extension	Myolite 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Morton's Ext.	U.skiv <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Rev. Morton's	U.skiv <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Dynamic Wedge	3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
	Sulcus Length	Sulcus Extension	Myolite 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Morton's Ext.	U.skiv <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Rev. Morton's	U.skiv <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Foot Cookie Ext.	Myolite 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
Partial Foot Toe Filler	Toe Filler <input checked="" type="checkbox"/>	5 th <input type="checkbox"/> 4 th <input type="checkbox"/> 3 rd <input type="checkbox"/> 2 nd <input type="checkbox"/> 1 st <input type="checkbox"/>		

E OFFLOADING PADS & CUSHIONS OPTIONAL

	Forefoot	Met Balance <input checked="" type="checkbox"/>	5 th <input type="checkbox"/> 4 th <input type="checkbox"/> 3 rd <input type="checkbox"/> 2 nd <input type="checkbox"/> 1 st <input type="checkbox"/>	
		Met Punch <input checked="" type="checkbox"/>	5 th <input type="checkbox"/> 4 th <input type="checkbox"/> 3 rd <input type="checkbox"/> 2 nd <input type="checkbox"/> 1 st <input type="checkbox"/>	
		Met Pad 2-4	4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Met Bar 1-5	4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Dancer's Pad	4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
	Midfoot	Cuboid Offload	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Rearfoot	Heel Cushion	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)
	Device Length Cushion	Heel Spur Pad	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Myolite Layer	4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
	Device Length Cushion	Plastazote Layer	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	

F FRAME MODIFICATIONS OPTIONAL

	Forefoot	1 st Ray Cut Out	Full <input type="checkbox"/> 65° <input type="checkbox"/> 45° <input type="checkbox"/>	
		5 th Ray Cut Out	65° <input type="checkbox"/> 45° <input type="checkbox"/>	
	Midfoot	Navicular B. Out	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		5 th Button Out	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
	Rearfoot	Heel Lift	6 <input type="checkbox"/> 4.5 <input type="checkbox"/> 3 <input type="checkbox"/> 1.5 <input type="checkbox"/> (mm)	
		Heel Lift Tapered	26 <input type="checkbox"/> 22 <input type="checkbox"/> 18 <input type="checkbox"/> 14 <input type="checkbox"/> 10 <input type="checkbox"/> (mm)	
		Heel Cup Depth	30 <input type="checkbox"/> 24 <input type="checkbox"/> 18 <input type="checkbox"/> 12 <input type="checkbox"/> 6 <input type="checkbox"/> (mm)	

Special Instructions:

