

# Drop Foot AFO: XTERN SIZE CHART



## AT-X TRIAL BRACE ASSESSMENT TOOL XTERN ADJUSTABLE TRIAL BRACE (AFO)

Men Shoe Size	Women Shoe Size	XTERN AFO & AT-X Size Selection	AT-X TRIAL BRACE ASSESSMENT TOOL XTERN ADJUSTABLE TRIAL BRACE (AFO)		
			AT -X	AT-X	AT-X
			LARGE	SMALL	PEDI
US-M	US-W	SIZE	AT-X FRONT LENGTH ADJUSTEMENT		
16		LARGE	61		
15.5		LARGE	57		
15		LARGE	53		
14.5		LARGE	50		
14		LARGE	46		
13.5		LARGE	42		
13		LARGE	38		
12.5		LARGE	34		
12		LARGE	31		
11.5		LARGE	28		
11		LARGE	25		
10.5	11.5	MEDIUM	22		
10	11	MEDIUM	19		
9.5	10.5	MEDIUM	16	52	
9	10	MEDIUM	13	49	
8.5	9.5	MEDIUM	11	46	
8 Y	8	MEDIUM	7	43	
7.5 Y	7.5	MEDIUM	3	40	
7 Y	7	8	0	37	
6.5 Y	6.5	7.5		34	
6 Y	6	7		30	
5.5 Y	5.5	6.5		26	
5 Y	5	6		22	
4.5 Y	4.5	5.5		18	36
4 Y	4	5		13	32
3.5 Y	3.5	4.5		9	29
3 Y	3	4		5	24
2.5 Y	2.5	3.5		1	21
2 Y					18
1.5 Y					16
1 Y					13
13 C					9
12.5 C					4
12 C					1
11 C					-3
10.5 C					---
10 C					---

Y: YOUTH SHOE SIZE (6-10 YEARS)  
C: CHILDREN SHOE SIZE (1-5 YEARS)

● CHOOSE SIZE FROM  
● PATIENT FOOT WIDTH



Call us: 1-877-767-3338  
Fax: 1-866-919-9268

www.kevinorthopedic.com  
hello@kevinorthopedic.com

Feel Change  
for a Better You

## CLINICIAN

1 REQUIRED

Account ..... Location .....

PO Number .....

Clinician .....

Clinician Email .....

## XTERN Drop Foot AFO Order form

Date: ..... / ..... / .....  
MM DD YYYY

Rush order due date:

..... / ..... / .....  
MM DD YYYY

Contact me to review Rx

**KEVIN**  
ORTHOPEDIC

Call us: 1-877-767-3338  
Fax: 1-866-919-9268  
www.kevinorthopedic.com  
hello@kevinorthopedic.com

## PATIENT

2 REQUIRED

Patient's Email .....

First Name .....

Last Name .....

DOB ..... / ..... / .....  
MM DD YYYY

Sex:  M  F Shoe Size

Weight ..... Height .....

Ship to Patient

Street Address .....

.....

City ..... State ..... Zip .....

## ORDER

3 REQUIRED

### SELECT XTERN AFO SIZE AND FRONT SHOE SECTION LENGTH

View Size Chart on reverse page for instructions on how to find both numbers

#### STEP 1

##### INDICATE XTERN AFO SIZE + QUANTITY

Based on shoe size. Use size chart on reverse page. The XTERN is not side specific. (If bilateral, indicate 2 below)

Small fits US M 2.5 - 8  
US W 4 - 9.5

Medium fits US M 7.5 - 10.5  
US M 9 - 12

Large fits US M 11 - 16  
US M 12.5 ▶



#### STEP 2


##### CHOOSE FRONT SHOE SECTION LENGTH

Based on size chart and clinical demo evaluation.

mm Small (0-55mm)

mm Medium (0-45mm)

mm Large (0-70mm)



#### REIMBURSEMENT

##### L1951

Ankle foot orthosis, spiral, (institute of rehabilitative medicine type), plastic or other material, prefabricated, includes fitting and adjustment.

##### L2820

Addition to lower extremity orthosis, soft interface for molded plastic, below knee section.

## OPTIONAL ACCESSORIES

### ANKLE STABILIZATION STRAP

Qty.



L2270 - Addition to lower extremity, varus/valgus correction ('t') strap, padded/lined or malleolus pad

### EXTENSION STOPPER KIT

Qty.



Visit the Kevin Orthopedic YouTube Channel for simple Patient Assessment guidelines ▶

Use QR App and focus the camera on your phone on this code to take you directly to video tutorial.



Special Instructions (including multiple product orders): .....