

# Offloading Pads & Cushions - *Midfoot*

Information Provided By Kevin Orthopedic Institute



## Cuboid Offloading Pad

Cuboid offloading pad of Myolite material

### FUNCTION:

- Offloads the cuboid

### CLINICAL INDICATIONS:

- Prominent cuboid bone
- Painful cuboid

A cuboid offload is a rectangular-shaped Myolite padding that is placed on the superior lateral surface of an orthotic's frame and sits beneath the distal calcaneus, proximal to the cuboid.

Cuboid Offload 6 4.5 3 1.5 (mm)



Inferior view  
(transparent)

6mm



4.5mm



3mm



1.5mm



Lateral view  
(transparent)

**Notes:** All illustrations and diagrams are of right foot  
Colors on illustrations are for visual purposes and will vary on final product