

1 CLINICIAN

Account Location

PO Number

Clinician

Clinician Email



Tel: 1-800-496-0987
 Fax: 1-866-919-9268
 www.kevinrootmedical.com
 hello@kevinrootmedical.com



Men's Shoes

[ORDER FORM PAGE 1 OF 2](#)

Date:/...../.....
MM DD YYYY

Contact me to review Order Form

Rush order due date:
/...../.....
MM DD YYYY

2 PATIENT

Patient's Email

First Name

Last Name

DOB/...../..... Sex: M F Weight Height Shoe Size

Ship to Patient

Street Address

City State Zip

3 MEN'S SHOES

SIZE 7 8 9 10 11 12 13 14

WIDTH Medium



Revere Durban Oiled Black

- Adjustable hook-and-loop straps.
- Leather uppers and linings.
- Removable comfort footbed with built-in arch support allows the insertion of custom orthotics if required.
- Matching strap extensions included to enhance fit.
- Slip-resistant outsole.
- Suitable for narrow widths.



Revere Montana 2 Oiled Black

- Adjustable hook-and-loop straps.
- Removable comfort footbed with build in arch support allows the insertion of custom orthotics if required.
- Matching strap extensions included to enhance fit.
- Slip-resistant outsole.



Revere Durban Whiskey

- Adjustable hook-and-loop straps.
- Leather uppers and linings.
- Removable comfort footbed with built-in arch support allows the insertion of custom orthotics if required.
- Matching strap extensions included to enhance fit.
- Slip-resistant outsole.
- Suitable for narrow widths.



Revere Montana 2 Whiskey

- Adjustable hook-and-loop straps.
- Removable comfort footbed with build in arch support allows the insertion of custom orthotics if required.
- Matching strap extensions included to enhance fit.
- Slip-resistant outsole.

4 FIT L3000 ORTHOTIC TO REVERE SANDAL


Select Modifications on reverse side or write special instructions

Order Custom Orthotic Fit to Sandals Bilateral Right Left

PATIENT

First Name

Last Name



Men's Shoes

Date:/...../.....
MM DD YYYY

[ORDER FORM PAGE 2 OF 2](#)

ANTI-PRONATION REACTION & POSTS - RIGHT **B** **LEFT - ANTI-PRONATION REACTION & POSTS**

	Varus Forefoot Posts	Balance FF to RF VR Extrinsic (Bar) VR Met Wedge	8° 6° 4° 2° 90° 6° 4° 2° 0° 6° 4° 2°			90° 2° 4° 6° 8° 0° 2° 4° 6° 2° 4° 6°	Balance FF to RF VR Extrinsic (Bar) VR Met Wedge	Varus Forefoot Posts	
	Midfoot	Raise Med Arch Scaphoid Pad VR Cuboid Pad Medial Flange	6 4.5 3 1.5 (mm) 6 4.5 3 1.5 (mm) 6 4.5 3 1.5 (mm) High Medium Low			(mm) 1.5 3 4.5 6 (mm) 1.5 3 4.5 6 (mm) 1.5 3 4.5 6 Low Medium High	Raise Med Arch Scaphoid Pad VR Cuboid Pad Medial Flange	Midfoot	
	Varus Rearfoot Posts	VR Intrinsic VR Extrinsic Heel Skive	6° 4° 2° 0° 6° 4° 2° 0° 6 4 2 (mm)			0° 2° 4° 6° 0° 2° 4° 6° (mm) 2 4 6	VR Intrinsic VR Extrinsic Heel Skive	Varus Rearfoot Posts	

ANTI-SUPINATION REACTION & POSTS - RIGHT **C** **LEFT - ANTI-SUPINATION REACTION & POSTS**

	Valgus Forefoot Posts	Balance FF to RF VG Extrinsic (Bar) VG Met Wedge	8° 6° 4° 2° 90° 6° 4° 2° 0° 6° 4° 2°			90° 2° 4° 6° 8° 0° 2° 4° 6° 2° 4° 6°	Balance FF to RF VG Extrinsic (Bar) VG Met Wedge	Valgus Forefoot Posts	
	Midfoot	Lower Med Arch Lateral Flange VG Frame Filler VG Cuboid Pad	6 4.5 3 1.5 (mm) High Medium Low 6° 4° 2° 0° 6 4.5 3 1.5 (mm)			(mm) 1.5 3 4.5 6 Low Medium High 0° 2° 4° 6° (mm) 1.5 3 4.5 6	Lower Med Arch Lateral Flange VG Frame Filler VG Cuboid Pad	Midfoot	
	Valgus Rearfoot Posts	VG Intrinsic VG Extrinsic	6° 4° 2° 0° 6° 4° 2° 0°			0° 2° 4° 6° 0° 2° 4° 6°	VG Intrinsic VG Extrinsic	Valgus Rearfoot Posts	

FOREFOOT EXTENSIONS - RIGHT **D** **LEFT - FOREFOOT EXTENSIONS**

	Toe Length	Toe Extension Morton's Ext. Rev. Morton's Dynamic Wedge	Myolite 4.5 3 1.5 (mm) U.skiv 4.5 3 1.5 (mm) Skived 4.5 3 1.5 (mm) 3 1.5 (mm)			(mm) 1.5 3 4.5 Myolite (mm) 1.5 3 4.5 U.skiv (mm) 1.5 3 4.5 Skived (mm) 1.5 3	Toe Extension Morton's Ext. Rev. Morton's Dynamic Wedge	Toe Length	
	Sulcus Length	Sulcus Extension Morton's Ext. Rev. Morton's Foot Cookie Ext.	Myolite 4.5 3 1.5 (mm) U.skiv 4.5 3 1.5 (mm) Skived 4.5 3 1.5 (mm) Myolite 4.5 3 1.5 (mm)			(mm) 1.5 3 4.5 Myolite (mm) 1.5 3 4.5 U.skiv (mm) 1.5 3 4.5 Skived (mm) 1.5 3 4.5 Myolite	Sulcus Extension Morton's Ext. Rev. Morton's Foot Cookie Ext.	Sulcus Length	
Partial Foot Toe Filler	Toe Filler	5 th 4 th 3 rd 2 nd 1 st			1 st 2 nd 3 rd 4 th 5 th	Toe Filler	Partial Foot Toe Filler		

OFFLOADING PADS & CUSHIONS - RIGHT **E** **LEFT - OFFLOADING PADS AND CUSHIONS**

	Forefoot	Met Balance Met Punch Met Pad 2-4 Met Bar 1-5 Dancer's Pad	5 th 4 th 3 rd 2 nd 1 st 5 th 4 th 3 rd 2 nd 1 st 4.5 3 1.5 (mm) 4.5 3 1.5 (mm) 4.5 3 1.5 (mm)			1 st 2 nd 3 rd 4 th 5 th 1 st 2 nd 3 rd 4 th 5 th (mm) 1.5 3 4.5 (mm) 1.5 3 4.5 (mm) 1.5 3 4.5	Met Balance Met Punch Met Pad 2-4 Met Bar 1-5 Dancer's Pad	Forefoot	
	Midfoot	Cuboid Offload	6 4.5 3 1.5 (mm)			(mm) 1.5 3 4.5 6	Cuboid Offload	Midfoot	
	Rearfoot	Heel Cushion Heel Spur Pad	6 4.5 3 1.5 (mm) 6 4.5 3 1.5 (mm)			(mm) 1.5 3 4.5 6 (mm) 1.5 3 4.5 6	Heel Cushion Heel Spur Pad	Rearfoot	
Device Length	Cushion	Myolite Layer Plastazote Layer	4.5 3 1.5 (mm) 6 4.5 3 1.5 (mm)			(mm) 1.5 3 4.5 (mm) 1.5 3 4.5 6	Myolite Layer Plastazote Layer	Device Length Cushion	

FRAME MODIFICATIONS - RIGHT **F** **LEFT - FRAME MODIFICATIONS**

	Forefoot	1 st Ray Cut Out 5 th Ray Cut Out	Full 65° 45° 65° 45°			45° 65° Full 45° 65°	1 st Ray Cut Out 5 th Ray Cut Out	Forefoot	
	Midfoot	Navicular B. Out 5 th Button Out	6 4.5 3 1.5 (mm) 6 4.5 3 1.5 (mm)			(mm) 1.5 3 4.5 6 (mm) 1.5 3 4.5 6	Navicular B. Out 5 th Button Out	Midfoot	
	Rearfoot	Heel Aperture Heel Lift Heel Lift Tapered	With Visco Plug No Plug 6 4.5 3 1.5 (mm) 26 22 18 14 10 (mm)			No Plug With Visco Plug (mm) 1.5 3 4.5 6 (mm) 10 14 18 22 26	Heel Aperture Heel Lift Heel Lift Tapered	Rearfoot	
Frame Attributes	Frame Filler	Rigid Semiflexible Myolite			Myolite Semiflexible Rigid	Frame Filler	Frame Attributes		

Special Instructions/Special Device: