

## Valisafe DentaCheck Protein Test (10 vials)

Order Details: **9VAL-PROTEINTEST**

The Denta-Check protein test is a rapid test that has the ability to detect residual proteins left behind on the surfaces of hard to clean dental and surgical instruments. An easy to read colour change gives you your final result within 10 seconds. Denta-Check is suitable for use with washer disinfectors, ultrasonics and even manual cleaning.

This fast, residual protein detection system helps you to identify that the load has been cleaned efficiently and can be used as a routine check. ISO15883-1 6.10.3 *Cleaning Efficacy Test 2* states a residual proteinaceous contamination is an **optional supplement to the mandatory daily 'visual' tests** to check for the presence of proteins, and during Performance Qualification (see Table A.1 1.3).


|                            |   |
|----------------------------|---|
| <b>Name:</b>               | Valisafe Denta-Check  |
| <b>Conformity:</b>         | BS EN ISO 15883, HTM 2030                                     |
| <b>Time</b>                | Result within 10 Seconds                                      |
| <b>Storage conditions:</b> | Refridgerated (4°C)   |
| <b>Lot no:</b>             | yes   |
| <b>Exp date:</b>           | yes   |
| <b>Correct symbols:</b>    | yes in accordance with the standard                           |
| <b>Indicator:</b>          | Diagnostic dye-binding vial                                   |
| <b>Expiry date:</b>        | 24 months refridgerated<br>6 months non-refridgerated         |
| <b>Packaging:</b>          | 10 tests per pack including swabs<br>+ 1 protein control test |

### Protein test which detects left behind residual on the Instruments

- Accurate results in 10 seconds
- No incubation required
- ISO15883 Compliant
- Clear colour change
- Easy to use
- No test more sensitive
- Cost effective
- Shelf life up to 2 years – when refrigerated
- Control test included to ensure batch has not been compromised



**Shelf Life** - 2 Years from date of manufacture if refrigerated, or 6 months if not refrigerated. Non refrigerated storage should not exceed 25deg C. Typically 12-18 months remaining at time of sale.



12 Juna Drive, Malaga, WA 6090      PO Box 1810, Malaga, WA 6944

T (08) 9244 4628

orders@sts-group.com.au

www.stshealth.com.au