



Ardo range of Breastpumps



Carum - Hospital grade



Calypso Double Plus



Single Electric Calypso



Amaryll Manual




Ardo Alyssa




Breast Shell gauge

A breast shell that does not fit correctly can prevent efficient expressing and cause pain and sore nipples. Please use the gauge to determine the correct size of breast shell for you*.


The Breast Shell size is correct if: the nipple moves freely in the funnel and can follow the rhythmic movements of the pump. No part, or very little, of the areola is in the funnel. The milk flows and the breast feels uniformly soft after expression.




For more details call 1300 GO ARDO (1300 46 2736) or visit www.ardo.com.au




22 mm




26 mm



28 mm



31 mm



36 mm

*A tolerance of 1 mm has been allowed. This gauge should only be used with the guidance of a healthcare professional.