

# PREvention<sup>®</sup>

EVERYDAY CONFIDENCE

Achieve **superior disinfection** of your Bellabaci massage cups with Prevention Disinfectants!

Our Disinfectant Wipes and Ready-to-Use Spray ensure safety, efficacy, and ease of use, empowering you with the confidence you deserve every day.



Active ingredient  
Hydrogen Peroxide ..... 0.5%  
Inert Ingredients ..... 99.5%  
Total ..... 100%

# BELLABACI CUPS CLEANING & DISINFECTION PROTOCOL

As a massage therapist, maintaining the cleanliness and hygiene of your massage cups is paramount to ensuring the well-being of your clients. Follow this step-by-step protocol to effectively clean and disinfect your Bellabaci cups, providing a safe and enjoyable massage experience.

**MATERIALS NEEDED:** Prevention Wipes or Ready-to-Use Spray, Soap or Mild Detergent, Warm Water, Soft Brushes or Cleaning Pads, Clean Towel, Disposable Gloves

## 1 Cleaning Steps

Put on disposable gloves to protect your hands.

Fill a basin or sink with warm water and add a small amount of soap or mild detergent.

Submerge the cups in the soapy water, ensuring all surfaces are covered.

Use soft cloth to gently scrub the exterior and interior of the cups.

Rinse the cups thoroughly with clean, warm water to remove any soap residue.

Wipe the cups with a clean, dry cloth.

The soapy water solution should be disposed of once cleaning is completed.

## 2 Disinfecting Steps

Gently wipe each cup evenly using Prevention Wipes or spray with Ready-to-Use Spray.

Allow the cups to remain wet for **1 minute** to kill bacteria, viruses and fungi.

Gently wipe the cups with a clean, damp cloth to remove any remaining detergent residue.

Now your disinfected massage cups are ready for your next client.

These instructions apply only to dry cupping equipment.

\*Refer to the product label or reference sheet for the full list of contact times, efficacy claims and use directions. **DISCLAIMER:** You must read and follow the instructions on the label for our product and use the product for the particular applications specified on the product label. To the full extent permitted under applicable law, Virox Technologies Inc. shall not be liable for any damages whatsoever, arising out of or otherwise in connection with any failure to follow those instructions or any use of its products which is not specified on the product's label. The instructions on the product label are not a warranty, express or implied, and Virox Technologies, Inc. expressly disclaims any implied warranties, including any warranty of fitness for a particular purpose. Under applicable law, some of the above waivers or disclaimers of warranty, or exclusions of liability for certain types of damages may not apply to you. In such cases, our liability will be limited to the fullest extent permitted by applicable law. Images are for illustrative purposes – please consult your state board for best practice guidelines and use of personal protection equipment. AHP® and Design and Accelerated Hydrogen Peroxide® and Design are trademarks of Diversey, Inc.