


Submitting a KLOwen Case through EasyRx

1. Log into your EasyRx Account, then click on "Create Rx."



The screenshot shows the EasyRx dashboard interface. At the top, there is a navigation bar with links for Home, Contact Support, Messages, and a toggle for Show tips. Below this is a secondary navigation bar with links for Dashboard, Patients, EasyRx 3D, Templates, Stats, Scripts, and Labs. The main content area features a sidebar on the left with various filters and a central table of prescriptions. A red box highlights the 'CREATE RX' button in the sidebar, and a red arrow points to it from the text above.

easyrx Orthodontics

dashboard needed by date awaiting feedback in-house lab aligner tracking

CREATE RX

Favorite Prescriptions

All Submitted (EasyRx Default View) SAVE VIEW AS SET AS DEFAULT PRINT VIEW

| Submitted Date | Last 90 Days | Office(All) | Doctor(All) | Tag(All) | Laboratory(All) | Received: All | Delivered: No | Search this queue | | | |
|----------------|--------------|-------------|-------------|----------|-----------------|---------------|---------------------------------|-------------------|-------|-------|------------------|
| ID | Patient | S | Tags | Dsc | Appt. Office | Doctor | Laboratory | Date Needed | Rec'd | Del'd | |
| 3595857 | John Doe | S | ● | 📅 | | Doctor | ULab Systems | 1/7/22 | N | N | view |
| 3595850 | Jane Doe | S | ● | 📅 | | Doctor | ULab Systems | 12/27/21 | N | N | view |
| 3593345 | | S | ● | 📅 | | | Ortho Orthodontics In-House Lab | 12/22/21 | N | N | view check-in |
| 3589358 | | CK | ● ● | 📅 | | | ULab Systems | 12/30/21 | N | N | view |
| 3588891 | | CK | ● ● | 📅 | | | Orthodont Lab - ODL | 12/28/21 | N | N | view |
| 3569496 | | CK | + | 📅 | | | KLOWen Braces | 1/11/22 | N | N | view |
| 3533181 | | CK | + | 📅 | | | KLOWen Braces | 1/4/22 | N | N | view |
| 3499129 | | CK | ● ● | 📅 | | | ULab Systems | 11/24/21 | N | N | view |
| 3489095 | | CK | ● ● | 📅 | | | Orthodont Lab - ODL | 12/2/21 | N | N | view |

Show 10

2. Select "KLOwen Braces" as the lab, and fill out patient info.

The screenshot shows the EasyRx Universal Lab Prescription Form. The form is divided into several sections:

- Lab:** A dropdown menu labeled "Select ..." is highlighted with a red box and a red arrow.
- Patient — Prescription Information:** This section contains several fields:
 - First Name** and **Last Name** (highlighted with a red box)
 - Date Needed** (highlighted with a red box)
 - Doctor**, **Office**, and **Ship To** (highlighted with a red box)
 - Scan/Imp. By** and **Submitted By** (highlighted with a purple box and a purple arrow)
- CREATE** button: A blue button located below the patient information fields.
- Standard Price:** A section showing a price of \$0.00.
- Parts Used:** A table listing various orthodontic parts and their corresponding arches.
- Comments:** A section for adding comments, including a "SAVE COMMENT" button.

A text overlay in the center of the form reads: "Choose a lab and enter a patient name to create a prescription."

*Note: The fields "Scan/Imp. By" and "Submitted By" applies to practices that have multiple users in one EasyRx Account for case tracking within the practice.

3. Add case-specific details in "Comments - Special Instructions."

Examples: extractions, missing teeth, other appliances in treatment, retain spaces for implants, supernumerary teeth, and anything different from a standard setup according to your typical preferences.

The screenshot displays the 'easyrx Universal Lab Prescription Form' interface. At the top, navigation links include Home, Contact Support, Messages, Hide, Show tips, and Orthodontics. A secondary navigation bar contains Dashboard, Patients, Easy Rx 3D, Templates, Stats, Scripts, Labs, and a search icon. The main form is divided into several sections:

- Lab:** A dropdown menu set to 'KLOWEN Braces' with the KLOWEN logo below it.
- Patient - Prescription Information:** Fields for Test (two instances), Month, Day, Year, Patient ID, Allergies, and Rx History. Appointment details include App't Date (12/24/2021), Date Needed (12/23/2021), Doctor, Scan/Imp. By, App't Time, Time Needed (optional), Office, Submitted By, Code, and Ship To.
- Rx #3597389:** Displayed on the right side.
- Library:** A section on the left with a search bar and a list of categories: Parts, Appliances, Templates, Favorites, Removables, Positioner, Sleep Therapy, Arch Development, Retention, Spring Aligner Appliances, Fixed Retention, Implants, Dentures, General Dental, Functionals, Fixed, Splints, Herbst, and Restorative.
- Standard Price:** \$0.00.
- Canvas:** Two dental arch diagrams labeled 'UPPER' and 'LOWER' with 'R' and 'L' markers.
- PARTS USED:** A tabbed interface with 'TEMPLATES' selected. It shows a list of parts with checkboxes for 'Upper' and 'Lower'.
- COMMENTS - SPECIAL INSTRUCTIONS:** A section at the bottom right, highlighted with a red box, containing a text area for 'Add a Comment', a 'SAVE COMMENT' button, and the 'EasyComments' logo.

4. Attach patient scans & photos in the "Upload" section.

- For iTero, enter "iTero Case Order ID" under the "Model Source" section, or click "Look Up and Insert iTero Case" to browse for patient scans (see next step).

The screenshot shows a web form for setting up an iTero case. At the top, there are tabs: PARTS USED, RESETS/MISSING, **TEMPLATES**, MODELS, ALIGNERS, IUD, and BANDS. Below these is a section for "Practice Assigned Templates". To the right is a "COMMENTS - SPECIAL INSTRUCTIONS" section with a text area, a "SAVE COMMENT" button, and the "EasyComments" logo. A red arrow points from this section down to the "UPLOAD" section. The main form is divided into three columns. The left column is "MISCELLANEOUS" with fields for "Outbound Tracking Number", "Coupon Code", "Consultation Call", and "OrthoAnalyzer P...". The middle column is "MODEL SOURCE" and is highlighted with a red box. It contains a radio button selected for "Use iTero Scan", a text field for "iTero Case Order ID", a "VALIDATE" button, and the "MyItero" logo. Below this is a "Look Up and Insert iTero Case" link and four other radio button options: "I am sending plaster models or impressions", "Use the digital model attached to this script", "Use the digital model that is in patient file", and "Other source" with a text field. The right column is "UPLOAD" and is also highlighted with a red box. It contains a text area with the instruction: "To attach files, drag and drop files here or select [files from your computer](#)". At the bottom of the form are two large buttons: "SAVE" (blue) and "SIGN & SUBMIT" (green). Below the buttons are links for "Terms and Conditions" and "Privacy", and a copyright notice: "Copyright 2009-2021 easyRx LLC".

iTero

SAVE

SIGN & SUBMIT

Terms and Conditions | Privacy

Copyright 2009-2021 easyRx LLC

- Minimum File Requirement = intra-oral scans + photo for patient verification.
- Other Preferred Files = all intra-oral photos, pano, cephs, CBCTs, etc.

5. When clicking "Look Up and Insert iTero Case", find the desired iTero scan and click "attach scan" to the right of it.

*This will attach all iTero files associated to that patient within MyAlignTech.

The screenshot shows the MyAlignTech interface with a modal window titled "iTero Scans". The modal contains a table of iTero scans and a search bar. A red arrow points to the "attach scan" link in the table.

iTero Scans

[Check for New Scans](#)

Search this queue


| Patient | Case Order ID | Status | Case Type Description | Scan Date | |
|---------|---------------|-----------|-----------------------|-----------|-----------------------------|
| | 77570573 | Completed | iRecord | 12/20/21 | attach scan |
| | 77569779 | Completed | iRecord | 12/20/21 | attach scan |
| | 77558211 | Completed | iRecord | 12/20/21 | attach scan |
| | 77546679 | Completed | iRecord | 12/20/21 | attach scan |
| | 77531620 | Completed | iRecord | 12/20/21 | attach scan |
| | 77530288 | Completed | iRecord | 12/20/21 | attach scan |
| | 77220744 | Completed | iRecord | 12/17/21 | attach scan |
| | 77215201 | Completed | iRecord | 12/16/21 | attach scan |
| | 77211504 | Completed | iRecord | 12/16/21 | attach scan |
| | 77208066 | Completed | iRecord | 12/16/21 | attach scan |

1 2 3 4 5 ... 326 NEXT

SAVE SIGN & SUBMIT

Terms and Conditions | Privacy
Copyright 2009-2021 easyRx LLC

6. Click "Sign & Submit" to finish.



The screenshot shows a web form with several sections. At the top, there are tabs: PARTS USED, RESETS/MODIFYING, **TEMPLATES**, MODELS, ALIGNERS, IUD, and BANDS. Below these tabs is a section titled "Practice Assigned Templates" which is currently empty. To the right of this section is a "COMMENTS - SPECIAL INSTRUCTIONS" box with a text area labeled "Add a Comment", a "SAVE COMMENT" button, and the "EasyComments" logo. Below the "Practice Assigned Templates" section is a "MISCELLANEOUS" section with input fields for "Outbound Tracking Number", "Coupon Code", "Consultation Call", and "OrthoAnalyzer P...". To the right of this is a "MODEL SOURCE" section with a radio button selected for "Use iTero Scan", a "VALIDATE" button, and a "MyItero" link. Below this are four other radio button options: "I am sending plaster models or impressions", "Use the digital model attached to this script", "Use the digital model that is in patient file", and "Other source". To the right of the "MODEL SOURCE" section is an "UPLOAD" section with a text area and a link to "files from your computer". At the bottom of the form are two buttons: a blue "SAVE" button and a green "SIGN & SUBMIT" button. A red arrow points from the top left towards the "SIGN & SUBMIT" button, which is also highlighted with a red rounded rectangle. At the very bottom of the page, there is a footer with the text "Terms and Conditions | Privacy" and "Copyright 2009-2021 easyRx LLC".

Restorative

PARTS USED | RESETS/MODIFYING | **TEMPLATES** | MODELS | ALIGNERS | IUD | BANDS

Practice Assigned Templates

COMMENTS - SPECIAL INSTRUCTIONS

Add a Comment

SAVE COMMENT

EasyComments

MISCELLANEOUS

Outbound Tracking Number

Coupon Code

Consultation Call

OrthoAnalyzer P...

MODEL SOURCE

☒ Use iTero Scan **VALIDATE** MyItero

Look Up and Insert iTero Case

☐ I am sending plaster models or impressions

☐ Use the digital model attached to this script

☐ Use the digital model that is in patient file

☐ Other source

UPLOAD

To attach files, drag and drop files here or select [files from your computer](#)

SAVE

SIGN & SUBMIT

Terms and Conditions | Privacy

Copyright 2009-2021 easyRx LLC