



Rehabilitation Collection


# Balto® Flexor



- Achilles tendon and ligament injuries
- Tarsal bone fractures
- Tarsus (hock) instability
- Hyperextension
- Hyperflexion
- Arthritis/ Arthrosis

The Balto® Flexor can be used as an alternative to surgery if your dog is not eligible for other options, such as age or previous injuries. The brace is designed to constrain joint movement while allowing an adequate range of motion (ROM), functioning as a rehabilitation support.

## Diagram



## Sizes

ITEM#	DISTANCE A-B		WEIGHT	
BTFLEXORXS-S	3.15 - 4"	8 cm - 10 cm	22 - 33 lbs	10-15 Kg
BTFLEXORM-L	4" - 5.5"	10 cm - 14 cm	33 - 99 lbs	15-45 Kg
BTFLEXORXL-XXL	5.5" - 7"	14 cm - 18 cm	99+ lbs	45+ Kg