



Rehabilitation Collection

# Balto<sup>®</sup> Hock

## Feline

- Achilles Tendon/Ligament Injury
- Tarsal Bone Fractures
- Tarsus (Hock) Instability
- Hyperextension or Hyperflexion
- Arthritis/Arthrosis



Addressing hock pain, lameness and neuromuscular problems, the Balto<sup>®</sup> Hock brace for cats also aides achilles injuries and orthopedic conditions, neurological conditions, pre/post-surgical immobilization, and contralateral support.

### Diagram



### Sizes

ITEM#	SIZE	WEIGHT	A	B	C
BTHOCKXXS	XXSW	5 - 10 lbs	3.75"	3"	1.75" - 2.75"
		2-5 kg	9.5 cm	7.5 cm	4.45 cm - 7 cm
BTHOCKXS	XS	10 - 20 lbs	4.25"	3.25"	2.25" - 3"
		5-8 kg	11 cm	8 cm	5.7 cm - 7.6 cm

Measurement A: Circumference above the hock  
 Measurement B: Circumference below the hock  
 Measurement C: Height of the brace