



Rehabilitation Collection


Balto® Joint

- Front leg trauma
- Minor carpus bone fractures
- Carpus joint instability
- Carpal inflammation
- Arthritis/ arthrosis



The Balto® Joint comfortably stabilizes the affected area, reducing pain and swelling while the animal heals. It can also provide protection for wraps and bandages from licking or biting. By removing the rigid splints, the brace can also be used as a compression brace for acute or chronic arthritis and arthrosis.

Diagram



Sizes

ITEM#	SIZE	WEIGHT	A	B
BTCARPALXS	X-Small	9 - 17.5 lbs	2"	2.75"
		4 - 8 kg	5 cm	7 cm
BTCARPALS	Small	17.5 - 55 lbs	3"	3"
		8 - 25 kg	7.6 cm	7.6 cm
BTCARPALM	Medium	55 - 110 lbs	4"	4.75"
		25 - 35 kg	10 cm	12 cm
BTCARPALL	Large	110 - 130 lbs	4.5"	5"
		35 - 50 kg	11.5 cm	12.7 cm
BTCARPALXL	X-Large	130+ lbs	5.25"	6.25"
		50+ kg	13 cm	16 cm

Measurement A: Circumference of metacarpus
 Measurement B: Height of the brace