



Rehabilitation Collection


# Balto® Pull



- Proprioceptive deficits
- Sciatic nerve damage
- Carpal flexural deformity/ Dorsiflexion
- Paw soreness
- Hyperflexion of phalanges
- Neurological conditions affecting paw placement
- Knuckling

The Balto® Pull addresses neurological disorders that cause severe proprioceptive deficits, such as dorsiflexion, or when the dog is unable to sense the position of the limb on the ground and tends to stand on the top of its paw, and damage to the sciatic nerve.

## Diagram



## Sizes

ITEM#	SIZE	WEIGHT	
BTPULLXXS	XX-Small	5 - 10 lbs	2 - 5 kg
BTPULLXS	X-Small	10 - 20 lbs	5 - 9 kg
BTPULLS	Small	20 - 45 lbs	9 - 20 kg
BTPULLM	Medium	45 - 85 lbs	20 - 39 kg
BTPULLL	Large	85+ lbs	39+ kg