



Rehabilitation Collection


Balto® Pull

- Proprioceptive deficits
- Sciatic nerve damage
- Carpal flexural deformity/ Dorsiflexion
- Paw soreness
- Hyperflexion of phalanges
- Neurological conditions affecting paw placement



The Balto® Pull addresses neurological disorders that cause severe proprioceptive deficits, such as dorsiflexion, or when the dog is unable to sense the position of the limb on the ground and tends to stand on the top of its paw, and damage to the sciatic nerve.

Diagram



Sizes

ITEM#	SIZE	WEIGHT
BTPULLXXS	XX-Small	5 - 10 lbs 2 - 5 kg
BTPULLXS	X-Small	10 - 20 lbs 5 - 9 kg
BTPULLS	Small	20 - 45 lbs 9 - 20 kg
BTPULLM	Medium	45 - 85 lbs 20 - 39 kg
BTPULLL	Large	85+ lbs 39+ kg

Measurement: Weight as indicated by a scale



Rehabilitation Collection

Balto® Pull


Feline



- Feline hyperflexion or hyperextension of phalanges
- Cases of “drop foot” or knuckling in cats
- Neurological conditions affected paw placement

The Balto® Pull for felines addresses neurological disorders that cause severe proprioceptive deficits in cats. This feline Achilles tendon brace is a great solution for cats suffering from limited mobility following injuries to the hock/Achilles tendon. This breathable anchor strap, made with non-neoprene material, provides stability and support to your cat’s rear legs without compromising natural movement.

Diagram



Sizes

ITEM#	SIZE	WEIGHT
BTPULLCATSXXXS	XX-Small	5 - 10 lbs. 2 - 5 kg
BTPULLCATSXS	X-Small	10 - 20 lbs. 5 - 9 kg

Measurement: Weight as indicated by a scale