



**Model: 58/3007**

**Anti-Slip Transfer Board with No Hand Holds**

Note: Persons with limited hand function can hook thumb into hand holes to assist board placement and removal.

We also recommend use of ADI's Push & Transfer gloves for those with limited grip and hand function.

**Features**

Sliding boards are used to promote independence and safety during transfers for those with limited upper body strength and impaired balance.

We also recommend our transfer boards as a transfer-assist tool for therapist / caregiver, to provide the needed assistance “to bridge the gap” for safer and more bio mechanically friendly transfer.

- Anti-Slip material on the bottom side keeps the board in place during transfers
- Tapered ends allow for easier placement
- Made from durable, lightweight hardwood
- Size: 74 x 20 x 1.2cm (29” x 8” x 0.5”)