



Model: 58/3006

Anti-Slip Transfer Board with Two Vertical Hand Holds

Note: Persons with limited hand function can hook thumb into hand holes to assist board placement and removal.

We also recommend use of ADI's Push & Transfer gloves for those with limited grip and hand function.

Features

Sliding boards are used to promote independence and safety during transfers for those with limited upper body strength and impaired balance.

We also recommend our transfer boards as a transfer-assist tool for therapist / caregiver, to provide the needed assistance "to bridge the gap" for safer and more bio mechanically friendly transfer.

- Anti-Slip material on the bottom side keeps the board in place during transfers
- Tapered ends allow for easier placement
- Made from durable, lightweight hardwood
- Size: 74 x 20 x 1.2cm (29" x 8" x 0.5")