

Oxy User's Manual



1- Power switch
Turn the device ON or OFF.

2- LCD Display
Displays the current test results.

3- Battery compartment
It works with two batteries AAA

Oxy User's Manual

How it works

1- Turn on the device

Press the power switch.

2- Connect to device

Open the SYHM app

2.1 - Android/ iOS apps:

- Open the app, select the language and then click OK.
- Click on Oximeter.
- Your device must be shown under Available Devices.
- Click on the device name.

2.2 - Windows / Mac OS X apps:

- Connect the USB bluetooth dongle to your pc
- Open the app and click on "Add Device +"
- Your device must now appear on the list.
- Click on the device name.

3- Insert the finger

Open the clamp insert a finger into the device.

4- Send the test results to the device

The application is now showing a notification indicating that the device is ready to receive the results.

Enter the desired value and click on "Send" button.

5- Check the test results

Test results appear on the LCD display.

Recharging the Battery

Remove the screws of the battery compartment.

Put the two batteries into battery compartment in correct polarities.