

SYSTEM 3.5 FOR SMALL FRAGMENTS




www.arzzt.com

- ▶ Distal Tibia Locking Plate
- ▶ Self-tapping Cortex Screw
- ▶ Stardrive Self-tapping Locking Screw
- ▶ Cancellous Screw



“Tapered end” for submuscular plate insertion, improving tissue viability.

- ▶ Manufactured in titanium Ti6Al4V-ELI
- ▶ Manufactured in Stainless Steel 316LS
- ▶ Extra-articular and simple intra-articular distal tibial fractures
- ▶ Distal tibial fractures, percutaneous or reducible by limited arthrotomy
- ▶ Distal tibial fracture extending into the diaphyseal area
- ▶ Limited contact profile. Design reduces plate to-bone contact, limiting vascular trauma
- ▶ Available for left and right tibias
- ▶ Shaft with 6 to 12 combination locking/compression holes
- ▶ 8 Distal locking holes accept $\varnothing 3.5\text{mm}$ locking
- ▶ Shaft $\varnothing 3.5$ locking, $\varnothing 3.5\text{mm}$ cortex, or $\varnothing 4.0\text{mm}$ cancellous bone screws in the compression portion
- ▶ Proximal and distal holes for $\varnothing 2.0$ mm Kirschner wires
- ▶ Distal tab for optional medial malleolus screw accepts $\varnothing 3.5\text{mm}$ locking, $\varnothing 2.7\text{mm}$ locking. The tab can bend, contour or cut if required.
- ▶ Anatomically contoured; plate is twisted and bent to fit distal Tibia



Approved by



System 3.5 for Small Fragments

Code 08070001

SYSTEM 3.5 FOR SMALL FRAGMENTS



www.arzzt.com

Code 08070001

Distal Tibia Locking Plate



Code	Side	Holes	Width
03141101	Right	6	138.73 mm
03141103	Right	8	164.73 mm
03141105	Right	10	191.75 mm
03141107	Right	12	214.74 mm
03141109	Right	14	244.75 mm
03141110	Left	6	138.73 mm
03141112	Left	8	164.73 mm
03141114	Left	10	191.75 mm
03141116	Left	12	214.74 mm
03141118	Left	14	244.75 mm

Self-tapping Cortex Screw



Ø 3.5 mm

Code	Length	Code	Length
02051201	10 mm	02051210	28 mm
02051202	12 mm	02051211	30 mm
02051203	14 mm	02051212	32 mm
02051204	16 mm	02051213	34 mm
02051205	18 mm	02051214	36 mm
02051206	20 mm	02051215	38 mm
02051207	22 mm	02051216	40 mm
02051208	24 mm	02051217	45 mm
02051209	26 mm	02051218	50 mm

Stardrive Self-tapping Locking Screw



Ø 2.7 mm

Ø 3.5 mm

Code	Length	Code	Length
02094203	14 mm	02091201	10 mm
02094204	16 mm	02091202	12 mm
02094205	18 mm	02091203	14 mm
02094206	20 mm	02091204	16 mm
02094207	22 mm	02091205	18 mm
02094208	24 mm	02091206	20 mm
02094217	45 mm	02091207	22 mm
		02091208	24 mm
		02091209	26 mm
		02091210	28 mm
		02091211	30 mm
		02091212	32 mm
		02091213	34 mm
		02091214	36 mm
		02091215	38 mm
		02091216	40 mm
		02091217	45 mm
		02091218	50 mm

Cancellous Screw



Ø 4.0 mm

Code	Length	Code	Length
02061201	10 mm	02061210	28 mm
02061202	12 mm	02061211	30 mm
02061203	14 mm	02061212	32 mm
02061204	16 mm	02061213	34 mm
02061205	18 mm	02061214	36 mm
02061206	20 mm	02061215	38 mm
02061207	22 mm	02061216	40 mm
02061208	24 mm	02061217	45 mm
02061209	26 mm	02061218	50 mm

Contact

202 Timberloch PI Suite 200
The Woodlands, TX, 77380
contact@arzzt.com