


ETCO₂ / O₂ NASAL SAMPLING CANNULA

DEHP & LATEX FREE | SINGLE PATIENT USE | SOFT PRONGS


- Designed to deliver oxygen supply and sampling of exhaled CO₂ gas simultaneously for more accurate EtCO₂ readings.
- Designed to prevent mixing of fresh oxygen with CO₂
- Crush-resistant tubing prevents occlusion and kinking
- Compatible with Microstream™, Capnostream™ and Oridion™ capnography monitors.
- Available in 2meters and 3meters length.



Designed for oxygen supply and sampling of exhaled CO₂ gas simultaneously.



SIZE	ADULT ONLY	CODE	PRICE
2METER	BOX/100	LB1430-2M	\$249
3METER	BOX/100	LB1430-3M	\$299



CONTACT US FOR MORE INFORMATION



Superior products, quality and service.
WWW.LSRHEALTHCARE.COM.AU

☎ 1300 859 767

✉ ENQUIRIES@LSRHEALTHCARE.COM.AU