

Khelix™ Fixed Curve Diagnostic Catheter

Features

- Flush electrode rings allow near-zero step height
- Cable platform allows multiple ring configurations

Tip: **Active platinum-iridium tip**

French Size: **6F**

Length: **115cm**

Product Code	Poles	Curve	Spacing
5SSF-ANN-4A252	4	Josephson	2-5-2 mm
5SSF-FNN-4A252	4	Cournand	2-5-2 mm
5SSF-DNN-4A252	4	Damato	2-5-2 mm
5SSF-YNN-4A555	4	HIS	5-5-5 mm

Product Code	Poles	Curve	Spacing
5SSF-ANN-6A252	6	Josephson	2-5-2 mm
5SSF-FNN-6A252	6	Cournand	2-5-2 mm
5SSF-DNN-6A252	6	Damato	2-5-2 mm
5SSF-YNN-6A555	6	HIS	5-5-5 mm



Khelix™ Decapolar Steerable Diagnostic Catheter

0°-180° Deflection

Features

- Flush electrode rings allow near-zero step height
- Cable platform allows multiple ring configurations
- Reliable push/pull mode of deflection

Type: **SUD**

Tip: **Active platinum-iridium tip**

Curves: **D/F/G**

French Size: **6F**

Length: **115cm**

Product Code

SS1-G0U-10A6222

SS1-G0U-10A6252

SS1-G0U-10A6282

SS1-G0U-10A6555


Spacing

2-2-2 mm

2-5-2 mm

2-8-2 mm

5-5-5 mm



Khelix™ DuoDecapolar Steerable Diagnostic Catheter

Features

- Flush electrode rings allow near-zero step height
- Superior signal quality demonstrated through clinical study¹
- Helical winding cable platforms allows multiple ring configurations
- Reliable push/pull mode of deflection

Type: **SUD**

Tip: **Passive**

Curves: **P/Q**

French Size: **6F**

Length: **115cm**

Product Code

SS2-P0U-20P6282

SS2-P0U-20P62102

SS2-P0U-20P62122

SS2-P0U-20P62132


Spacing

2-8-2 mm

2-10-2 mm

2-12-2 mm


2-13-2 mm



0°-270° Deflection



Khelix™ Variable Loop Diagnostic Catheter



Features

- Flush electrode rings allowing near-zero step height
- Loop safety release mechanism
- Integrated torque tube
- Cable platform allows multiple ring configurations
- Reliable push/pull mode of deflection
- Compatible with 8.0 - 8.5F steerable and non-steerable long introducers

Tip: **Inactive non-bulbous dome**

Diameter: **25-15mm**

Loop: **4.5F**

Curve: **B**

French Size: **7F**

Length: **115cm**

Product Code

Poles

Spacing

SC2-B2515U-10P7888

10

8-8-8 mm


SC2-B2515U-20P7262

20

2-6-2 mm

Khelix™ Fixed Loop Steerable Diagnostic Catheter

0°-180° Deflection



Features

- Flush electrode rings allowing near-zero step height
- Loop safety release mechanism
- Integrated torque tube
- Cable platform allows multiple ring configurations
- Reliable push/pull mode of deflection
- Compatible with 8.0 - 8.5F steerable and non-steerable long introducers



Tip: **Inactive non-bulbous dome**

Poles: **10**

Loop: **4.5F**

Curve: **B**

French Size: **7F**

Length: **115cm**

Product Code

Diameter

Spacing

SC1-B15U-10P7555

15mm

5-5-5 mm

SC1-B20U-10P7555

20mm

6-6-6 mm

SC1-B25U-10P7888

25mm

8-8-8 mm

HERES WHY YOU SHOULD BUY THE
KHELIX

○ AUSTRALIAN MADE AND OWNED

○ BETTER PRICING

○ REPORTED IMPROVED SIGNAL ACQUISITION

