

VOIR BOUGIE®

Endotracheal tube introducer with Safety Markers




5mm x 70cm

FOR USE WITH ENDOTRACHEAL TUBES 6.0MM ID AND LARGER IN ADULT PATIENTS


COLORED SAFETY MARKINGS | LATEX FREE | SINGLE USE

The Patented safety markings are a unique feature of the disposable Voir Bougie.

Adroit Surgical's introducers are manufactured from low density polyethylene which is both firm and flexible for easy insertion.



The Voir Bougie has brightly marked colored bands at both ends that allows rapid accurate insertion to the correct depth by the user and helps to keep the bougie in the same position during use to avoid patient injury or inadvertent removal of the bougie.



LSR HEALTHCARE PTY LTD
WWW.LSRHEALTHCARE.COM.AU
ENQUIRIES@LSRHEALTHCARE.COM.AU

