

SWEDE-O® WALKING BOOTS

- Helps immobilize, support, and protect the ankle and foot
- Rocker bottom promotes natural gait
- Cushioned sole for shock absorption
- Wide footbed and deluxe liner provides a comfortable fit
- Lightweight yet durable construction
- Uprights are anatomically contoured for a better fit
- Bilateral design fits right or left foot



STUDY REFERENCE

¹ <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5721918/>

Available From:



Short

Tall

HELPS SUPPORT ANKLE & FOOT

Swede-O® Walking Boots help immobilize the ankle and is ideal for severe ankle sprains, stable ankle fractures, post surgery, and post cast recovery. Studies¹ show that bearing weight on a stable injury helps facilitate recovery, allowing the patient to return to normal activity at a faster rate.

Designed with a non slip rocker bottom and low profile heel to promote natural gait. Shock absorbing, closed cell foam insole and heel cushion absorb and dissipate the impact of heel strike when walking. A wide footbed provides increased stability and accommodates for edema and bandage while a deluxe foam liner provides comfort.

Swede-O Walking Boot is durable, lightweight, and includes fitting pads to accommodate extra space around the ankle, heel, or leg. Anatomically contoured uprights provide a comfortable fit.

SPECIFICATIONS

Product Selection: Swede-O® Walking Boot, Black

Short: UTL-1130 | Tall: UTL-1131

Product Weight:

Short: 2.3 lbs | Tall: 2.8 lbs

Short Dimensions: Small: 10.5" x 11" | Medium: 11.5" x 11.5" | Large: 13" x 13"

Tall Dimensions: Small: 15.5" x 10.5" | Medium: 17.5" x 11.5" | Large: 18.5" x 13"

Materials: pvc vinyl, polyethylene

Case Quantity: 6

Case Dimensions:

Short: 17" x 14" x 18" | Tall: 23" x 14" x 18"

Case Weight:

Short: 15 lbs | Tall: 18 lbs

Latex: This product is not made with natural rubber latex

Country of Origin: Made in China

California Residents Only:

⚠️ WARNING: Cancer and Reproductive Harm – www.P65Warnings.ca.gov