

## SWEDE-O® PEDIATRIC WALKING BOOT

- Helps immobilize, support, and protect the ankle and foot
- Rocker bottom promotes natural qait
- Cushioned sole for shock absorption
- Wide footbed and deluxe liner provides a comfortable fit
- Lightweight yet durable construction
- Bilateral design fits right or left foot





STUDY REFERENCE

https://www.ncbi.nlm.nih.gov/pmc/ articles/PMC5721918/

Available From:



## **HELPS SUPPORT ANKLE & FOOT**

Swede-O® Pediatric Walking Boot helps immobilize the ankle and is ideal for severe ankle sprains, stable ankle fractures, post surgery, and post cast recovery. Studies<sup>1</sup> show that bearing weight on a stable injury helps facilitate recovery, allowing the patient to return to normal activity at a faster rate.

Designed with a non-slip rocker bottom and low profile heel to promote natural gait. Shock absorbing, closed cell foam insole and heel cushion absorb and dissipate the impact of heel strike when walking. A wide footbed provides increased stability and accommodates for edema and bandage while a deluxe foam liner provides comfort.

Swede-O Pediatric Walking Boot is durable, lightweight, and includes fitting pads to accommodate extra space around the ankle, heel, or leg.

## **SPECIFICATIONS**

**Product Selection: Swede-O® Pediatric** Walking Boot, Gray UTL-1132

Product Weight: .8 lbs

**Short Dimensions:** 

Small: 7.5" x 6" Medium: 8.5" x 7" Large: 10" x 8"

XLarge: 12" x 9"

Materials: pvc vinyl, polyethylene

Case Quantity: 3

Case Dimensions: 12" x 8" x 10"

Case Weight: 4 lbs

Latex: This product is not made with natural rubber latex

Country of Origin: Made in Vietnam

California Residents Only:

▲WARNING: Cancer and Reproductive Harm — www.P65Warnings.ca.gov