

SOFT COMFORT CORPAK™ HOT OR COLD THERAPY PACKS

- Soft surface quickly releases hot or cold therapy
- Remains flexible, comfortably contours to painful areas
- Frost free exterior is gentle on the skin
- Store in freezer, heat in microwave
- Nontoxic, biodegradable fill maintains a therapeutic temperature for 20 minutes
- Various sizes can accommodate any treatment area
- Customized packs can be used as a promotional item (call for details)

Tri-Sectional



Medium



Small



Cervical



Large



X-Large



Available From:



HELP RELIEVE PAIN

Soft Comfort CorPaks™ Hot or Cold Therapy Packs are the perfect solution to help relieve headaches, ease pain, and help reduce inflammation associated with strains, sprains, and other soft tissue injuries. Conveniently store packs in the freezer for cold therapy the moment you need it! These packs stay soft and flexible, comfortably contouring to painful areas.

COLD THERAPY

Use the Soft Comfort CorPak cold to help reduce inflammation and swelling and provide soothing relief from muscle and joint pain, back, shoulder, neck, and knee discomfort. Use for sports injuries, post-workout recovery, and rehabilitation. Also ideal for fever reduction, chronic pain, bursitis, tendinitis, and arthritis.

HOT THERAPY

Use Soft Comfort Corpak hot to accelerate recovery and healing by improving circulation and blood flow to areas of discomfort. Increasing the temperature of the afflicted area even slightly can soothe and relax muscles, promote healing of damaged tissue, and increase muscle flexibility. These packs are also ideal for menstrual cramp relief when applied to the lower abdomen or lower back.

ORGANIC, NONTOXIC & BIODEGRADABLE

The pack's exterior is made of a proprietary blend of synthetic fibers and organic plant materials. Our Soft Comfort CorPaks contain more nontoxic, biodegradable fill than traditional packs to retain cold and heat longer, while helping reduce environmental impact. The durable RF welded seam provides a secure seal to help prevent leaks.

SPECIFICATIONS

Product Selection: Soft Comfort CorPak's™ Hot or Cold Therapy Packs

Materials: Food grade propylene glycol USP, Bentonite clay, Water
Cover: Wood pulp, Polyester, Polyolefin alloy film

Latex: This product is not made with natural rubber latex

Country of Origin: Made in USA

California Residents Only:
⚠️WARNING: Cancer and Reproductive Harm – www.P65Warnings.ca.gov

SOFT COMFORT CORPAK™ HOT OR COLD THERAPY PACKS

INSTRUCTIONS FOR USE

For cold therapy, place pack in freezer for a minimum of 1 hour prior to use.

For heat therapy, simply heat the pack in a microwave according to the heating instructions printed on the pack.

The pack will maintain therapeutic temperature for 20 minutes. When not in use, the pack may be stored in the freezer for convenience.

All Soft Comfort CorPak's are:



New!

Pictured Above:
X-Large Soft Comfort CorPak™ (12.5" x 18.5")

Available From:



Small

ACC-549
Product Dimensions: 3" x 5"
Product Weight: .25 lbs
Case Quantity: 12
Case Weight: 4 lbs
Case Dimensions: 11.5" x 11.5" x 3"



Medium

ACC-550
Product Dimensions: 6" x 10"
Product Weight: 1.4 lbs
Case Quantity: 16
Case Weight: 25 lbs
Case Dimensions: 12" x 12" x 8"



Large

ACC-551
Product Dimensions: 10" x 13"
Product Weight: 3.2 lbs
Case Quantity: 6
Case Weight: 21 lbs
Case Dimensions: 12" x 12" x 6"



Cervical

ACC-552
Product Dimensions: 6" x 20"
Product Weight: 1.88 lbs
Case Quantity: 10
Case Weight: 21 lbs
Case Dimensions: 12" x 12" x 6"



Tri-Sectional

ACC-554
Product Dimensions: 11" x 15"
Product Weight: 3.75 lbs
Case Quantity: 6
Case Weight: 26 lbs
Case Dimensions: 12" x 12" x 6"



New!

X-Large

ACC-555
Product Dimensions: 12.5" x 18.5"
Product Weight: 5.8 lbs
Case Quantity: 4
Case Weight: 25 lbs
Case Dimensions: 12" x 12" x 6"