


## ALGERBRUSH INSTRUMENTS ARE AVAILABLE IN RUMEX RANGE

### COMPLETE SET OF A POWER HANDLE, CHUCK, AND BURR

- LIGHTWEIGHT
- POCKET-SIZED
- BATTERY POWERED
- REPLACEMENT BURRS AVAILABLE



USED WORLDWIDE SINCE 1970  
CE MARKED SINCE 1998

### PTERYGIUM REMOVER

Diamond burr provides excellent smoothing of corneal surface after surgical removal of the tissue.

- 16-051-2.5** AlgerBrush with a 2.50 mm Diamond Round Fine Grit Burr
- 16-050-3.5** AlgerBrush with a 3.50 mm Diamond Round Medium Grit Burr
- 16-050-5.0** AlgerBrush with a 5.00 mm Diamond Disk-Shaped Medium Grit Burr



### RUST RINGS REMOVER

Used to remove foreign bodies and rust rings from patient's eye. Tungsten carbide coating ensures efficient corneal surface smoothing.

- 16-140** AlgerBrush Rust Ring Remover with a 1.00 mm Burr
- 16-141** AlgerBrush Rust Ring Remover with a 0.50 mm Burr



### REPLACEMENT BURRS

Diamond burrs for PTERYGIUM REMOVER  
Chuck and burr, non-sterile

- 16-051-2.5B** Round Fine Grit Burr for AlgerBrush, 2.50 mm
- 16-051-3.5B** Round Medium Grit Burr for AlgerBrush, 3.50 mm
- 16-052-5.0B** Disk-Shaped Medium Grit Burr for AlgerBrush, 5.00 mm

Tungsten burrs for RUST RING REMOVER  
Burr only, 5 per box, non-sterile

- 16-143B** Burr for AlgerBrush, 0.50 mm
- 16-142B** Burr for AlgerBrush, 1.00 mm

**Not for sale in Canada**