

EzAim™ SoftBite™

Adhesive **Cushioned** Bite Tab Holders




Comfort in Every Bite

Cushioned Tab
for Extra Patient Comfort

Shortened Design
for Greater Accessibility

Peel-and-Stick Adhesive



Compare to Schick AimRight™ and Dentsply FASTab™

The EzAim SoftBite sensor holder with peel-and-stick adhesive technology firmly adheres to digital sensors while allowing for extra patient comfort with a cushioned bite tab. The cushioned bite tab allows patients to comfortably bite harder on the tab for more precise imaging while having a shortened design for greater accessibility. Single-use, disposable bite tabs with no extra parts or setup needed.

- Shortened design for greater accessibility
- Cushioned tab for extra patient comfort
- No extra parts needed for set up
- Easily adjustable to any patient anatomy
- Strong adhesive allows for reliable positioning

ITEM#	DESCRIPTION	QUANTITY
131115	EzAim SoftBite Disposable Adhesive Cushioned Bite Tab Holders, Red	100/box

BUY 3, GET 1 FREE!




Scan to Request FREE Samples

© 2024 Pac-Dent, Inc. All rights reserved.

Pac-Dent®
Passion for Excellence














EzAim™

X-Ray Positioning System



The EzAim X-Ray Positioning System features disposable, self-adhesive sensor holders with peel-and-stick adhesive technology that firmly secures digital sensors to holders. The holders are color coded for easy identification and fits all types of sensors on the market. EzAim holders can be used with existing aiming rings and positioning arms. Patient-friendly design eliminates extra dimensions for setup. Can be managed with hemostats. Available with matching, autoclavable positioning arms and aiming rings.


BUY 3, GET 1 of the Same FREE!

ITEM#	DESCRIPTION	COMPARED TO	QUANTITY
Sensor Holder Tabs			
131109	Blue Anterior Sensor Holder Tab	Xios Plus 6176510	100/pack 
131110	Yellow Posterior Sensor Holder Tab	Xios Plus 6176528	100/pack 
131111	Red Bitewing Sensor Holder Tab	Xios Plus 6176536	100/pack 
131112	Green Endo Sensor Holder Tab	Xios Plus 6176544	100/pack 
131113	Red Bite Tab Sensor Holder	Schick AimRight B1073051	50/pack 
Positioning Arms			
131101	Blue Anterior Positioning Arm	XCP/BAI 540857	1/pack 
131102	Red Bitewing Positioning Arm	XCP/BAI 540927	1/pack 
131103	Yellow Posterior Positioning Arm	XCP/BAI 540858	1/pack 
131104	Green Endo Positioning Arm	XCP/BAI 550597	1/pack 
Aiming Rings			
131105	Blue Anterior Aiming Ring	XCP/BAI 540865	1/pack 
131106	Red Bitewing Aiming Ring	XCP/BAI 540934	1/pack 
131107	Yellow Posterior Aiming Ring	XCP/BAI 540860	1/pack 
131108	Green Endo Aiming Ring	XCP/BAI 550598	1/pack 
Intro Kit			
131114	4 x Positioning Arms (Anterior Blue, Bitewing Red, Posterior Yellow, and Endo Green), 4 x Aiming Rings (Anterior Blue, Bitewing Red, Posterior Yellow, and Endo Green), 125 Disposable Adhesive Sensor Holders (25 of each Anterior Blue, Bitewing Red, Posterior Yellow, Endo Green, and Bite Tab Red)		

Aiming Rings

Autoclavable


- Easy color coded identification system
- System helps to eliminate consuming reprocessing time



Positioning Arms

Autoclavable


- Compatible with existing aiming rings
- Autoclavable steel components



Sensor Holder Tabs

Disposable

- Peel-and-stick adhesive
- Universal digital sensor compatibility



① Anterior ② Bitewing ③ Posterior ④ Endo ⑤ Bite Tab



Scan to Request FREE Samples

© 2024 Pac-Dent, Inc. All rights reserved.

Pac-Dent®
Passion for Excellence