

Technical specs (tips)

iVac 35 tip

- Color: **green**
- O.D.: **0.35mm (± 0.02)**
- I.D.: **0.15mm (± 0.02)**
- Taper (D_0 - D_{16}): **0.025 (± 0.05) [$D_0=35$; $D_{16}=75$]**
- Total length: **30mm**
- Working length: **27mm**
- Depth marks at 18, 19, 20, and 23mm

iVac 50 tip

- Color: **yellow**
- O.D.: **0.50mm (± 0.02)**
- I.D.: **0.30mm (± 0.02)**
- Taper (D_0 - D_{16}): **0.025 (± 0.05) [$D_0=50$; $D_{16}=90$]**
- Total length: **30mm**
- Working length: **27mm**
- Depth marks at 18, 19, 20, and 23mm

Angled capillary tip

- O.D.: **0.60mm (± 0.05)**
- I.D.: **0.30mm (± 0.05)**
- Working length (after angle): **21mm**