


PAC-DENT PROPHY ANGLE FAMILY

Disposable Propy Angle Comparison Chart



Prophy Angle	AntiSplatr™	Mini 90™	Mighty Mini 90™	Mini Ergo™	ProAngle® EZ	ProAngle® Ergo	FUNimals™	ProAngle® Plus	ProAngle® Plus with Brush Cup	ProAngle® Plus Pedo	ProAngle®
MSRP	\$0.71/Each	\$0.74/Each	\$0.70/Each	\$0.74/Each	\$0.69/Each	\$0.59/Each	\$0.62/Each	\$0.41/Each	\$0.76/Each	\$0.38/Each	\$0.36/Each
Body Size	Mini	Mini	Mini	Mini	Extra Long	Regular	Regular	Regular	Regular	Mini	Regular
Body Type	90 Degree	90 Degree	90 Degree	Contra Angle	90 Degree	Contra Angle	90 Degree	90 Degree	90 Degree	90 Degree	90 Degree
Description	<ul style="list-style-type: none"> Effectively reduce splatter by up to 100% Patented wiper removes paste before clusters form and are flung during tooth polishing Golf ball dimples keep saliva attached to the exterior surface of the cup 	<ul style="list-style-type: none"> Compact head and slimmer profile allow superior visibility and handling Improved flare for better subgingival access Patented beveled gear design eliminates vibration and ensures smooth, quiet operation Pre-sterilized with Gamma radiation to reduce the risk of cross-contamination 	<ul style="list-style-type: none"> Compact, powerful head for increased visibility and improved access Improved flare for better subgingival access Patented beveled gear design delivers quiet and smooth power for enhanced patient comfort 	<ul style="list-style-type: none"> Ergonomically designed with the optimal amount of head tilt to reduce hand and wrist fatigue 20% shorter head and 25% slimmer neck profile for increased visibility and improved access Patented beveled gear design delivers quiet and smooth power for enhanced patient comfort 	<ul style="list-style-type: none"> Reduce the weight of the handpiece, which gives dentists and hygienists a much easier time on their wrists and hands Eliminates the costs of cleaning and disinfection for the nose cone With its patented beveled gear design, the angle runs smoothly with no overheating 	<ul style="list-style-type: none"> Ergonomic design reduces hand and arm The gearing structure of ProAngle Ergo makes it much quieter than other contra-angle DPAs Torque Cup minimizes splatter, improves efficiency by 15%, and offers better interproximal cleaning 	<ul style="list-style-type: none"> Built on the critically acclaimed ProAngle Plus, the FUNimals Disposable Propy Angles are the newest addition to the line of high quality Pac-Dent DPAs These FUN angles will reduce your young patients' anxiety while offering optimal performance 	<ul style="list-style-type: none"> Reinforced beveled gears to ensure only one angle is required for each procedure. Eliminates vibration for less noise, chatter, and heat Guaranteed no cup dropping 	<ul style="list-style-type: none"> Features the same mechanical design and technology of ProAngle Plus but with flat or tapered brush cups 	<ul style="list-style-type: none"> Features the same patented ProAngle design but with petite, latex-free ribbed cups Ideal for pediatric and small mouth patients and better accessibility to hard-to-reach areas Cups are made from rubber material to ensure they last the entire operation 	<ul style="list-style-type: none"> Patented prophy angle that employs the same design as ProAngle Plus to ensure each prophy procedure only requires one angle to finish Utilizes the same beveled gear design as ProAngle Plus for assured quality and performance
Configuration	100/pk, 500/pk	144/pk, 500/pk	144/pk, 500/pk	144/pk, 500/pk	100/pk, 500/pk	144/pk, 500/pk	100/pk, 500/pk	144/pk, 500/pk	100/pk	144/pk, 500/pk	144/pk, 500/pk
Available Colors	Pink, Blue	Turquoise	Green	Champagne	Blue	Teal, Blue, Lavender	Assorted animal prints	Purple, Teal, Blue	Purple, Teal	Purple, Teal	Purple
Compatible Handpieces	All	All	All	All	Pac-Dent ProMate™ E-Type Motor, Midwest® Shorty®/Rhino®, Star® Titan®	All	All	All	All	All	All
Cup Options	Soft, Firm	Soft, Firm, Pedo	Soft	Soft, Firm, Pedo	Soft, Firm, Tapered Brush	Super-Soft, Soft, Firm, Torque	Soft	Super-Soft, Soft, Firm, Torque	Flat, Tapered Brush Cup	Soft(Pedo), Firm(Pedo)	Soft, Firm
Ergonomically design	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Optimal Flare	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	N/A	Yes	Yes
Sterile	No	Yes	No	Yes	No	No	No	No	No	No	No
Latex-Free	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes



Scan QR code to request free samples