

# Guide to Face Shields Selection and Use

A variety of designs help meet the individual needs of healthcare workers



Choose the appropriate facial protectant device to meet your individual needs and anticipated level of exposure and risk.

## Visor & Face Shields

Pictured: iVisor™ Flip Visor + Shield Kit  
Also available in almond and pink

### iVisor™ Flip Visor + Shield Kit

Reusable lightweight visor and shield kit with innovative easy-flip face shield. Perfect for healthcare workers who require frequent access to his or her face. Features tinted overhead minishield for glare reduction and cinch cord for fast on-the-go adjustment.



## Face Shields

Pictured: iShield™ Pro Flip Face Shield

### iShield™ Pro Flip Disposable Face Shields

Disposable ultralight face shield with convenient flip-up barrier shield. Perfect for first responders who require the lightest full-coverage face protection. Features fog-free shield, vented sponge head cushion, and adjustable elastic headband.



## Visor & Face Shields

Pictured: iVisor™ Loupe Visor + Shield Kit  
Available in multiple colors and sizes

### iVisor™ Loupe Visor + Shield Kit

Ultra-comfortable reusable visor and shield kit designed with extra large breathing room and wide air vents for maximum breathability and airflow. Easily fits in front of optical accessories such as loupes, lights, and glasses.



## Face Shields

Pictured: iShield™ Pro Face Shield

### iShield™ Pro Disposable Face Shields

Disposable ultralight face shield designed to offer the highest user comfort. Features extra thick vented sponge cushion that absorbs sweat and provides excellent ventilation. Recommended for prolonged wear.



## Visor & Face Shields

Pictured: iVisor™ Visor + Shield Kit  
Available in multiple colors

### iVisor™ Visor + Shield Kit

Classic reusable shield and visor kit for full-face protection against sprays, splashes, and splatters of body fluids and hazardous materials. Clamp-style fit with fog-free shields and thick foam head bumper.



## Face Shields

Pictured: iShield™ Face Shield

### iShield™ Disposable Face Shields

Classic disposable ultralight face shield to provide wearers full frontal face protection from hazardous fluids. Vented sponge cushion and elastic headband ensure even pressure and eliminate stress points.



## Face Shield Drapes

Pictured: Armor™ Drape-It™ Face Shield Drapes

### Armor™ Drape-It™ Disposable Face Shield Drapes

Apply Drape-It to any facial protectant device to instantly increase coverage and protection.



## Ordering Information

	ITEM#	DESCRIPTION	COLOR	QTY
<b>iVisor Visor &amp; Shield Kit:</b> 3 face shields, 1 mini shield, 1 LC filter, 1 user manual				
<b>iVisor</b>	<b>V013B</b>	iVisor visor & shield kit with blue visor	Blue	1 kit/pk
	<b>V013PK</b>	iVisor visor & shield kit with pink visor	Pink	1 kit/pk
	<b>V013V</b>	iVisor visor & shield kit with violet visor	Violet	1 kit/pk
	<b>V013W</b>	iVisor visor & shield kit with white visor	White	1 kit/pk
	<b>V013T</b>	iVisor visor & shield kit with teal visor	Teal	1 kit/pk
	<b>V013BK</b>	iVisor visor & shield kit with black visor	Black	1 kit/pk
<b>iVisor Loupe Visor &amp; Shield Kit:</b> 3 face shields, 2 mini shields, 1 LC filter, 1 cinch cord, 1 user manual				
<b>iVisor Loupe</b>	<b>V201P-S</b>	iVisor Loupe visor & shield kit with pink visor, Small 8 cm	Pink	1 kit/pk
	<b>V201PW-S</b>	iVisor Loupe visor & shield kit with pearl white visor, Small 8 cm	Pearl White	1 kit/pk
	<b>V201B-M</b>	iVisor Loupe visor & shield kit with blue visor, Standard 9 cm	Blue	1 kit/pk
	<b>V201T-M</b>	iVisor Loupe visor & shield kit with teal visor, Standard 9 cm	Teal	1 kit/pk
	<b>V201P-M</b>	iVisor Loupe visor & shield kit with pink visor, Standard 9 cm	Pink	1 kit/pk
	<b>V201SB-M</b>	iVisor Loupe visor & shield kit with space blue visor, Standard 9 cm	Space Blue	1 kit/pk
	<b>V201PW-M</b>	iVisor Loupe visor & shield kit with pearl white visor, Standard 9 cm	Pearl White	1 kit/pk
	<b>V201SBK-M</b>	iVisor Loupe visor & shield kit with satin black visor, Standard 9 cm	Satin Black	1 kit/pk
	<b>V201T-L</b>	iVisor Loupe visor & shield kit with teal visor, Large 10 cm	Teal	1 kit/pk
	<b>V201SB-L</b>	iVisor Loupe visor & shield kit with space blue visor, Large 10 cm	Space Blue	1 kit/pk
<b>iVisor Flip Visor &amp; Shield Kit:</b> 3 face shields, 2 mini shields, 1 LC filter, 1 cinch cord, 1 user manual				
<b>iVisor Flip</b>	<b>V401T</b>	iVisor Flip visor & shield kit with teal flip visor	Teal	1 kit/pk
	<b>V401A</b>	iVisor Flip visor & shield kit with almond flip visor	Almond	1 kit/pk
	<b>V401P</b>	iVisor Flip visor & shield kit with pink flip visor	Pink	1 kit/pk
<b>Refill</b>	<b>V014</b>	6 x mini shield replacements, <i>iVisor use only</i>	Clear	6/pk
	<b>V015</b>	3 x disposable face shield replacements, <i>for iVisor, iVisor Loupe &amp; iVisor Plus</i>	Clear	3/pk
	<b>V016</b>	12 x disposable face shield replacements, <i>for iVisor, iVisor Loupe &amp; iVisor Plus</i>	Clear	12/pk
	<b>V017</b>	6 x small mini shield replacements, <i>for iVisor Loupe &amp; iVisor Plus</i>	Tinted	6/pk
	<b>V018</b>	6 x large mini shield replacements, <i>for iVisor Loupe &amp; iVisor Plus</i>	Tinted	6/pk
	<b>V019</b>	6 x light-cure filter replacements, <i>for iVisor, iVisor Loupe &amp; iVisor Flip</i>	Clear	6/pk
	<b>V020</b>	6 x mini shield replacements, <i>iVisor Flip use only</i>	Clear	6/pk
	<b>V021</b>	12 x disposable face shield replacements, <i>iVisor Flip use only</i>	Tinted	12/pk
	<b>V022</b>	3 x disposable face shield with soft edge guard, <i>for iVisor, iVisor Loupe &amp; iVisor Plus</i>	Clear	3/pk
	<b>V023</b>	3 x disposable face shield with soft edge guard, <i>for iVisor, iVisor Loupe &amp; iVisor Plus</i>	Tinted	3/pk
<b>iShield</b>	<b>S024</b>	Disposable face shields. Vented	Clear	24/pk
	<b>S026</b>	Disposable face shields. Vented	Clear	38/pk
<b>iShield Pro</b>	<b>S030</b>	Disposable face shields. Vented	Clear	38/pk
<b>iShield Pro Flip</b>	<b>S032</b>	Disposable face shields. Vented	Clear	38/pk
<b>Armor Drape-It</b>	<b>ADI-144</b>	Disposable face shield drapes	Clear	144/pk

# Face Shields Comparison

	iVisor	iVisor Loupe	iVisor Flip	iShield	iShield Pro	iShield Pro Flip
<b>Full-face Coverage</b>	Yes	Yes	Yes	Yes	Yes	Yes
<b>Reusable</b>	Yes	Yes	Yes	No	No	No
<b>Cold-sterilizable</b>	Yes	Yes	Yes	No	No	No
<b>Anti-fog Coating</b>	Yes	Yes	Yes	Yes	Yes	Yes
<b>Overhead Protection</b>	Yes	Yes	Yes	No	No	No
<b>Easy-flip</b>	No	No	Yes	No	No	Yes
<b>Weight</b>	Lightweight	Lightweight	Lightweight	Ultralight	Ultralight	Ultralight
<b>Latex-free</b>	Yes	Yes	Yes	Yes	Yes	Yes
<b>Comfort Fit</b>	Clamp style	Adjustable cinch cord	Adjustable cinch cord	Elastic headband	Elastic headband	Elastic headband
<b>Light-cure Filter</b>	Yes	Yes	Yes	No	No	No

Follow CDC guidelines for using personal protection equipment (PPE). Face shields act as a physical barrier to block sprays, splashes, and splatters of infectious fluids before they reach the face. However, due to the inherent design and shape of face shields this protection is not absolute as infectious particles and droplets may travel under the face shield and be breathed in by the wearer. Therefore, it is strongly recommended that face shields be worn in conjunction with tight-fitting face masks and protective eyewear to cover all the mucous membranes on one's face that are most susceptible to infection.

Pac-Dent, Inc. makes no expressed or implied guarantee that these products will fully protect the user from exposure to hazardous fluids or prevent the contracting of infectious diseases. It is the user's due diligence to select the appropriate personal protection equipment to meet the expected level of exposure in accordance with OSHA regulations and CDC guidelines.

2020 SUPER SPONSOR OF



ISO 13485:2016  
ISO 9001:2000

