

Guide to Protective Eyewear Selection and Use

A variety of designs help meet the individual needs of healthcare workers



Choose the appropriate protective eyewear to meet your individual needs and anticipated level of exposure and risk.

Disposable Eyewear

Pictured: iWear™ Disposable Eyewear

iWear™ Disposable Eyewear

Ultralight disposable protective eyewear for ocular protection against hazardous fluids. Features fog-resistant lenses that provide superior optical clarity in wet or steamy environments. Available with a snap-on or loop-through frame.




Disposable Eyewear

Pictured: iWear™ Ultra Eyewear

iWear™ Ultra Disposable Eyewear

Sleek disposable eye shield for the modern healthcare worker. Ultralight frame and ergonomic design allows iWear Ultra to fit comfortably over glasses. Features fog-resistant lenses and foldable legs.



Disposable Eyewear

Pictured: iWear™ Plus Eyewear

iWear™ Plus Disposable Eyewear

Ultralight disposable protective eyewear specially engineered for maximum user comfort during prolonged wear. Included silicone nose pad and neck strap offer exceptional comfort and convenience. Features fog-free lenses and anti-slip grooved legs.



Reusable Eyewear

Pictured: iWear™ Safety Glasses

iWear™ Safety Glasses

Reusable lightweight safety glasses offer stylish protection in any workplace environment. Feature UV-blocking lenses and flexible legs that diffuse temple pressure made from impact-resistant polycarbonate. Can be decontaminated via cold sterilization.



Reusable Goggles

Pictured: iWear™ Full Coverage Protective Goggles

iWear™ Full Coverage Protective Goggles

Lightweight full-coverage goggles protect eyes from sprays, splashes, dust, and flying debris. Adjustable elastic headband, soft nose piece, and comfort-seal design allow a tight yet comfortable seal on the face. Easily fits over glasses. Can be decontaminated via cold sterilization.



Ordering Information

	ITEM#	DESCRIPTION	COLOR	QTY
iWear Disposable Eyewear: 10 frames with assorted colors, 10 clear lenses, 1 user manual				
iWear	G1010	iWear Disposable Eyewear	Assorted	10/pk
	G1060	iWear Disposable Eyewear	Assorted	10/pk
iWear Ultra Disposable Eyewear: 10 frames with blue colors, 10 clear lenses, 1 user manual				
iWear Ultra	G2010	iWear Ultra Disposable Eyewear	Blue	10/pk
iWear Plus Disposable Eyewear: 10 frames with assorted colors, 20 clear lenses, 1 neck strap, 1 nosecone, 1 user manual				
iWear Plus	G1020	iWear Plus Disposable Eyewear	Assorted	10/pk
	G100	10 clear lenses refill	Clear	10/pk
	G005	50 frames with assorted colors & 50 clear lenses	Assorted	50/pk
	G006	10 bonding lens for curing	Orange	10/pk
Safety Glasses	SG012	iWear Safety Glasses	White	12/pk
	SG012BK	iWear Safety Glasses	Black	12/pk
Goggles	CX02	iWear Full Coverage Goggles	White	10/pk

Eyewear Comparison

	iWear	iWear Ultra	iWear Plus	Safety Glasses	Protective Goggles
Full Ocular Coverage	Yes	Yes	Yes	Yes	Yes
Reusable	No	No	No	Yes	Yes
Cold-sterilizable	No	No	No	Yes	Yes
Anti-fog Coating	Yes	Yes	Yes	Yes	Yes
Fits Over Glasses	Yes	Yes	Yes	No	Yes
Comfort-seal Fit	No	No	No	No	Yes
Weight	Ultralight	Ultralight	Ultralight	Lightweight	Lightweight
Silicone Nose Piece	No	Yes	No	No	Yes
Impact-resistant	Yes	Yes	Yes	Yes	Yes
Latex Free	Yes	Yes	Yes	Yes	Yes

Follow CDC guidelines for using personal protection equipment (PPE). Protective eyewear help shield the eyes against exposure to infectious droplets and airborne debris generated during patient procedures. However, this protection only extends to the ocular area. It is strongly recommended to wear additional PPE such as face shields and face masks to meet the expected level of exposure to pathogens and other hazardous materials.

Pac-Dent, Inc. makes no expressed or implied guarantee that these products will fully protect the user from exposure to hazardous fluids/particles or prevent the contracting of infectious diseases. It is the user's due diligence to select the appropriate personal protection equipment to meet the expected level of exposure in accordance with OSHA regulations and CDC guidelines.

2020 SUPER SPONSOR OF



ISO 13485:2016
ISO 9001:2000

