

RUSH FOOT[®] EVAQ8 COLLECTION

Joel Goltry, BK, RUSH FOOT HiPro EVAQ8


**SIMPLISTIC.
EFFECTIVE.
EFFICIENT.**

Finally, an unprecedented elevated vacuum system that seamlessly integrates effective vacuum with the same rugged simplicity of all the **RUSH Foot[®]** designs. The **RUSH Foot EVAQ8 Collection** is the easy solution to providing maximum vacuum efficiency in just a few quick steps. No batteries or complicated, mechanical pump required.

Unlike the competition, upwards of 22inHg can be achieved with this surprisingly simple, integrated design that stabilizes volume, increases linkage, and heightens proprioception through elevated vacuum suspension. The **RUSH Foot EVAQ8 foot Collection** is all-terrain, fully-submersible (fresh or salt water), free of bulky components, and requires minimal maintenance. The **RUSH FOOT EVAQ8 Collection** is available in the **RAMPAGE, HiPro, ROGUE 2, and RAMPAGE LP** models providing the perfect option for the majority of individuals seeking elevated vacuum suspension.

**Go anywhere. Do anything.
The RUSH Foot EVAQ8
Collection elevates your life.**

RUSH ROGUE[®] 2 EVAQ8 PROFILE INFORMATION

| Weight Limit | 365lbs./165kg | Heel Height | 3/8"/10mm | Size | Build Height** |
|-----------------|---------------|--|-----------------------|-----------|----------------|
| Size | 22cm - 29cm | Warranty | 36 Months | 22 - 24cm | 6.125" / 155mm |
| Foot Shell | Light or Dark | ***Codes | L5987, L5984 L5781 | 25 - 27cm | 6.625" / 168mm |
| Product Weight* | 1039g |  | | 28 - 29cm | 6.75" / 171mm |

RUSH HiPro[®] EVAQ8 PROFILE INFORMATION

| Weight Limit | 365lbs./165kg | Heel Height | 3/8"/10mm | Size | Build Height** |
|-----------------|---------------|------------------|--------------|-----------|----------------|
| Size | 22cm - 29cm | Warranty | 36 Months | 22 - 24cm | 6.5" / 165mm |
| Foot Shell | Light or Dark | †Suggested Codes | L5981, L5781 | 25 - 27cm | 6.875" / 175mm |
| Product Weight* | 930g | | | 28 - 29cm | 7.125" / 181mm |

RUSH RAMPAGE LP EVAQ8 PROFILE INFORMATION

| Weight Limit | 365lbs./165kg | Heel Height | 3/8"/10mm | Size | Build Height** |
|-----------------|---------------|------------------|-----------------------|-----------|----------------|
| Size | 22cm - 29cm | Warranty | 36 Months | 22 - 24cm | 4.5" / 114mm |
| Foot Shell | Light or Dark | †Suggested Codes | L5981, L5986 L5781 | 25 - 27cm | 4.75" / 121mm |
| Product Weight* | 815g | | | 28 - 29cm | 4.75" / 121mm |

RUSH RAMPAGE EVAQ8 PROFILE INFORMATION

| Weight Limit | 365lbs./165kg | Heel Height | 3/8"/10mm | Size | Build Height** |
|-----------------|---------------|------------------|-----------------------|-----------|----------------|
| Size | 22cm - 29cm | Warranty | 36 Months | 22 - 24cm | 6.625" / 168mm |
| Foot Shell | Light or Dark | †Suggested Codes | L5987, L5986 L5781 | 25 - 27cm | 7.15" / 181mm |
| Product Weight* | 860g | | | 28 - 29cm | 7.25" / 184mm |

CATEGORY CHART

| Weight | lbs | 0 - 105 | 106 - 140 | 141 - 175 | 176 - 210 | 211 - 245 | 246 - 280 | 281 - 315 | 316 - 350 | 351 - 365 |
|-------------------|--------|---------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| kg | 0 - 48 | 49 - 64 | 65 - 79 | 80 - 95 | 96 - 111 | 112 - 127 | 128 - 143 | 144 - 159 | 160 - 166 | |
| Low Activity | | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 |
| Moderate Activity | | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 |
| High Activity | | 1 (2) | 2 (3) | 3 (4) | 4 (5) | 5 (6) | 6 (7) | 7 (8) | 8 (9) | 9 (-) |
| High Impact | | 2 (3) | 3 (4) | 4 (5) | 5 (6) | 6 (7) | 7 (8) | 8 (9) | 9 (-) | - (-) |

† The listing of this code should not be considered a guarantee for coverage or a payment.

* Product weight based on size 26cm, Cat 4 foot module with spectra sock and foot shell.

**Build height based on size 23cm, 26cm, or 29cm Cat 4 foot module with spectra sock, foot shell, and 10mm heel height.

***PDAC Verified

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability. (RAMPAGE LP)

Call us toll-free at 855.450.7300 or visit us at proteorusa.com

f t YouTube in p i #HumanFirst

